



## ▶ Bottle Caps and Reservoirs

- Positive containment of laboratory chemicals and vapors
- Chemically inert wetted parts
- Robust and easy to use
- Anti-twist design
- GL45, GL32 and 38-430 sizes
- 2, 3 and 4 port caps

Safe, simple and economic liquid containment. Open solvent bottles can allow harmful solvent vapor to escape and chemicals to be spilled. Omnifit bottle caps and accessories help prevent these problems and give you the security of a safe and dependable solvent delivery system.

### Products:

- **T Series - Luer ported straight through tubing**
- **Q Series: 1/4"-28 UNF ported with optional valves**
- **Reservoirs**
- **Accessories for Bottle Caps and Reservoirs**



 **Bottle Caps & Reservoirs**

**HRoMalytic** +61(0)3 9762 2034 09  
AUSTRALIAN Distributors **ECHnology** [www.chromtech.net.au](http://www.chromtech.net.au)

[Home](#) | [Contact Us](#)

Copyright © 2009 Bio-Chem Fluidics.