

## Series CGV-130

Fixed needle: .028" x .005" x 3.25", bevel, open end  
or .028" x .008" x 3.25", bevel, open end

- ▶ **Side vent** – for easy cleaning and filling
- ▶ **Reinforced plunger** – eliminates plunger bending
- ▶ **Plunger guide in barrel** – smaller front section of plunger is confined within guide section
- ▶ **Leak-tight to 250 psi** – liquids and gases



#139023

**Add a side vent for  
easy cleaning and filling**

The CGV-130 features the same rugged construction as the CG-130, with the addition of a side vent that can be used for filling or cleaning the syringe. A heavy, reinforced plunger within an enlarged plunger guide section of the barrel virtually eliminates plunger bending. The smaller front section of the plunger is confined within the guide section even when the plunger is retracted, while the sturdy upper section absorbs the force of a rapid downstroke.

	10 µl	25 µl	50 µl	100 µl
With .005" ID needle:	139022	139023	139024	139025
With .008" ID needle:			139124	139125