

# Autosampler and OEM Syringes

## Autosampler Syringes

### Hewlett Packard

This 10 µl syringe has a sturdy plunger to stand up to continuous automatic operation. The PTFE plunger tip prevents blow-back or freeze-up.



	OD	ID	Length	Type	Gauge	P/N
HP 7670	.028" x	.005" x	2.0"	Bevel, open end	22	136022
HP 7671	.028" x	.005" x	2.0"	Bevel, open end	22	136022
HP 7673A	.019" x	.006" x	1.75"	Blunt, open end	26	101986
HP 7673A	.024" x	.005" x	1.75"	Blunt, open end	23	101987
HP 7673A	.019" x	.005" x	2.0"	Bevel, open end	26	101988

### Varian

We presently supply the syringes for the Varian GC 511, PTA30, PTA30WS, and Archon instruments.



	OD	ID	Length	Type	Gauge	P/N
511-5 µl	.024" x	.006" x	2"	Taper, blunt, & open end	23	102100
511-50 µl	.024" x	.006" x	2"	Taper, blunt, & open end	23	102102
511-1.5 µl	.028" x	.008" x	2"	Taper, blunt, & open end	22	102103
PTA				Luer fitting		030038
PTA30WS				Luer fitting		030038
Dynatrap				Luer fitting		030038
Archon				26 ml barrel with decal		502968