

Liquid Syringes

Features

Precision Sampling liquid syringes have all the standard features on page 2, with special features delineated in the chart below.

	Pressure rating	Removable needle	Fixed needle	Reinforced plunger	Plunger stop	Side Vent	Page number
Series C-160:	250 psi	✓	✓				14
Series C-120:	250 psi		✓				15
Series CG-130:	250 psi	✓	✓	✓			16
Series CGV-130:	250 psi		✓	✓		✓	17
Series D-140:	30 psi	✓		✓	✓		18
Mini-injector:			✓	✓	✓		19
Purge and Trap:		luer					22

Sizes Available

	Series C-160	Series C-120	Series CG-130	Series CGV-130	Series D-140	Mini-injector	Purge & trap
Page no:	14	15	16	17	18	19	22
.01 µl:						✓	
.10 µl:						✓	
1 µl:					✓		
2 µl:					✓		
5 µl:	✓	✓	✓		✓		
10 µl:	✓	✓	✓	✓	✓		
25 µl:	✓	✓		✓			
50 µl:	✓	✓		✓			
100 µl:	✓	✓		✓			
5 ml:							✓
10 ml:							✓