



VICI® Valco Instruments Co. Inc.
VICI AG International

HOME | PRODUCTS | CUSTOMER SERVICE | SEARCH

Syringes

PRODUCTS

- [What's new](#)
- [Actuators](#)
- [Calibration gas standards](#)
- [Columns](#)
- [Fittings](#)
- [Flow, pressure, and on/off controllers](#)
- [Gas purifiers](#)
- [GC accessories](#)
- [GC detectors](#)
- [HayeSep porous polymers](#)
- [Injectors/valves](#)
- [Liquid handling](#)
- [Portable gas monitors](#)
- [Syringes](#)
- [Tools](#)
- [Trace gas analyzers](#)
- [Tubing](#)

CUSTOMER SERVICE

- [Support](#)
- [Catalog request](#)
- [Ordering information](#)
- [Returns/repairs](#)
- [Credit reference form](#)
- [Contact us](#)
- [Privacy statement](#)
- [Terms of sale](#)

VICI

- [About us](#)
- [Photo tour](#)
- [Upcoming trade shows](#)

Replacement needles for gas syringes

from VICI Precision Sampling

Needles are manufactured to exacting tolerances from high precision stainless steel tube stock. Slip-on needles for Series A-2, C, and D syringes feature a needle extension through the hub that fits inside the PTFE tip on the syringe. This design eliminates the dead volume inherent in conventional hubbed needles.

The point styles and lengths listed below are standard, but we are willing to supply others. Let us know your requirements by [emailing](#) VICI Precision Sampling or calling toll-free (800) 828-1653.



1. Blunt tip, open end



2. Bevel, open end



3. Side port, taper

Product numbers

Needles are sold in packages of three.

Series	Size	Type	Product No.
Pressure-Lok needles			
A	.029" x .012" x 2.25"	2	913050
	.029" x .012" x 2.25"	3	913052
A-2	.029" x .005" x 2"	2	943050
	.029" x .012" x 2"	2	943051
	.029" x .012" x 2"	3	943052
C	.028" x .005" x 2"	2	943050
	.029" x .012" x 2"	2	943051
	.029" x .012" x 2"	3	943052
D	.029" x .012" x 2"	2	943051
	.029" x .012" x 2"	3	943052

Magnum	.035" x .023" x 2"	2	953050
	.035" x .023" x 2"	3	953052
Luer needles			
A-2	.028" x .006" x 2"	2	943060
	.028" x .012" x 2"	2	943061
	.028" x .012" x 2"	3	943062
C	.028" x .006" x 2"	2	943060
	.028" x .012" x 2"	2	943061
	.028" x .012" x 2"	3	943062
D	.028" x .012" x 2"	2	943061
	.028" x .012" x 2"	3	943062
Magnum	.035" x .023" x 2"	2	953061
	.035" x .023" x 2"	3	953062

Other replacement parts

TOP ▲

[Home](#) | [Products](#) | [Customer service](#) | [Contact](#) | [About VICI](#) | [Site map](#) | [Site index](#)

VICI[®], Valco[®], and Cheminert[®] are registered trademarks of Valco Instruments Co. Inc. and VICI AG
© 2008 VICI