



20167



24785

Standard Microliter Syringes

for Agilent 7673, 7683, 7693A, and 6850 Autosamplers

- Hamilton and SGE syringes are designed and tested to meet critical autosampler performance needs.
- Needle point styles are designed to withstand multiple, rapid injections through a septum.

Hamilton

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek cat.#
					Model	cat.#	
0.5 µL	ASRN	23s	1.71"	Agilent	7000.5	86276	ea. 24895
0.5 µL	ASRN	26s	1.71"	Agilent	7000.5	86274	ea. 21221
5 µL	ASRN	23s	1.71"	Agilent	75	87957	ea. 20172
5 µL	ASN	23s	1.71"	Agilent	75	87987	ea. 20168
5 µL	ASN	23s	1.71"	Agilent	75	87990	6-pk. 20170
5 µL	ASN	26s	1.71"	Agilent	75	87988	ea. 24592
5 µL	ASN	26s	1.71"	Agilent	75	87989	6-pk. 21230
5 µL	ASRN	23s-26s	1.71"	Agilent	75	87959	ea. 21223
5 µL	ASN	23s-26s	1.71"	Agilent	75	87993	ea. 24593
5 µL	ASN	23s-26s	1.71"	Agilent	75	87994	6-pk. 24594
10 µL	ASRN	23s	1.71"	Agilent	701	80357	ea. 20171
10 µL	ASN	23s	1.71"	Agilent	701	80387	ea. 20167
10 µL	ASN	23s	1.71"	Agilent	701	80390	6-pk. 20169
10 µL	ASRN	26s	1.71"	Agilent	701	80358	ea. 24597
10 µL	ASN	26s	1.71"	Agilent	701	80388	ea. 24595
10 µL	ASN	26s	1.71"	Agilent	701	80389	6-pk. 24599
10 µL	ASRN	23s-26s	1.71"	Agilent	701	80359	ea. 24598
10 µL	ASN	23s-26s	1.71"	Agilent	701	80393	ea. 24596
10 µL	ASN	23s-26s	1.71"	Agilent	701	80391	6-pk. 24600

SGE

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek cat.#
					Model	cat.#	
0.5 µL	R	23	42 mm	Cone	0.5BR-HP-0.63	000410	ea. 24791
0.5 µL	R	26	42 mm	Cone	0.5BR-HP-0.47	000400	ea. 24790
5 µL	F	23	42 mm	Cone	5F-HP-0.63	001810	ea. 24781
5 µL	F	23	42 mm	Cone	SK-5F-HP-0.63	001814	6-pk. 24783
5 µL	F	26	42 mm	Cone	5F-HP-0.47	001800	ea. 24780
5 µL	F	26	42 mm	Cone	SK-5F-HP-0.47	001804	6-pk. 24782
5 µL	F	23-26s	42 mm	Cone	5F-HP-0.63/0.47	001821	ea. 21210
5 µL	F	23-26s	42 mm	Cone	SK-5F-HP-0.63/0.47	001822	6-pk. 21214
10 µL	R	23	42 mm	Cone	10R-HP-0.63	002815	ea. 24795
10 µL	F	23	42 mm	Cone	10F-HP-0.63	002810	ea. 24785
10 µL	F	23	42 mm	Cone	SK-10F-HP-0.63	002814	6-pk. 24787
10 µL	R	26	42 mm	Cone	10R-HP-0.47	002805	ea. 24794
10 µL	F	26	42 mm	Cone	10F-HP-0.47	002800	ea. 24784
10 µL	F	26	42 mm	Cone	SK-10F-HP-0.47	002804	6-pk. 24786
10 µL	F	23-26	42 mm	Cone	10F-HP-0.63/0.47	002821	ea. 21212
10 µL	F	23-26	42 mm	Cone	SK-10F-HP-0.63/0.47	002822	6-pk. 21215
Large Volume							
50 µL	R	23	42 mm	Cone	50R-HP-0.63	004665	ea. 24799