

Useful Gadgets



Now available for
GL-45 bottles!



Waste Overflow Indicator

- Avoid messy pooling around mobile phase waste containers.
- Audible alarm instantly alerts user, preventing overflow.
- Compact, battery operated unit.
- Available for 4-liter and GL-45 solvent bottles.

The Restek Waste Overflow Indicator will help keep your mobile phase waste where it belongs—in the waste container! Compact, battery operated unit fits securely on solvent bottles and accommodates two waste streams. An audible alarm is given as the solvent waste container approaches capacity, giving you time to empty or change the container. Another innovative design from Restek!

Description

| | qty. | cat.# | price |
|--|-------|-------|-------|
| Waste Overflow Indicator for HPLC Systems, 4 Liter | ea. | 26543 | |
| Waste Overflow Indicator for HPLC Systems, GL-45 | ea. | 26550 | |
| Replacement AA Battery for the Waste Overflow Indicator | ea. | 26544 | |
| Replacement AA Batteries for the Waste Overflow Indicator | 3-pk. | 26545 | |

**No more burned
fingers!**



26262

The Claw and The Claw Holder Kit

- Easily removes hot liners from injection ports.
- 4mL vials (not included) can be replaced when dirty.

Never again will you have to burn your fingers removing a hot injection port liner. The Claw safely and cleanly removes liners, O-rings, or other small objects from the injection port. You can then place the hot objects in a clean 4 mL vial situated in the Claw holder until ready for reuse.



Description

| Description | qty. | cat. # | price |
|--|---------|--------|-------|
| The Claw | ea. | 26261 | |
| The Claw Holder Kit (includes The Claw and holder) | kit | 26262 | |
| WISP 48 Snap Seal Vial | 100-pk. | 24658 | |

Description
Sapphire Scribe


qty. **cat.#** **price**
ea. 20182



Shortix® Capillary GC Column Cutter

- Consistently make precise, clean, square cuts with a diamond blade.
- Built-in magnifier to verify square cut.
- Use with 0.25 mm ID to 0.53 mm ID tubing (0.78 mm OD maximum).

| Description | qty. | cat.# | price |
|---|-------------|--------------|--------------|
| Shortix Capillary GC Column Cutter | ea. | 23026 | |
| Maintenance Kit for Shortix Capillary GC Column Cutter (Includes: diamond cutting wheel, O-rings, and a tool to open the column cutter) | kit | 23027 | |





Stainless Steel Jet Reamers

A great tool for cleaning detector jets or removing stuck ferrules from other small orifices. Serrated design is optimal for removing silica deposits and other contaminants.

| Description | qty. | cat.# | price |
|-----------------------------|-------------|-------|-------|
| Stainless Steel Jet Reamers | 6-piece set | 20113 | |



Ferrule Removal Kit

The tapered tools in this kit have teeth designed to grip and remove fused silica adaptor ferrules that have become stuck in the fitting detail. Each kit has two tools: one for removing 1/32-inch adaptor ferrules and one for removing 1/16-inch adaptor ferrules.

| Description | Valco # | qty. | cat.# | price |
|---------------------|---------|------|-------|-------|
| Ferrule Removal Kit | FRK1 | kit | 20146 | |



Mini Hand Drill Set

Drill ferrules to the proper ID in seconds! Includes three drills, for use with 0.25, 0.32, and 0.53 mm ID (0.4, 0.5, and 0.8 mm OD) capillary columns.

| Description | qty. | cat.# | price |
|---------------------|-------------|-------|-------|
| Mini Hand Drill Set | 3-piece set | 20122 | |

1/16-Inch Tubing Cutter

- Produces square, smooth cuts in 1/16-inch tubing.
- Eliminates tubing distortion.
- Replaceable cutting wheel.



Description

1/16" Tubing Cutter

Replacement Cutting Wheels

qty.

ea.

3-pk.

cat.#

20192

20185

price

Ridgid® Tubing Cutter

- Excellent for cutting 1/8- or 1/4-inch metal tubing.*
- Compact size is ideal for tight spaces.
- Replaceable cutting wheel.



Description

Ridgid Tubing Cutter for 1/8" or 1/4" metal tubing

Replacement Cutting Wheels

qty.

ea.

2-pk.

cat.#

23011

23012

price

1/16-Inch Tubing Cutting Pliers

- Ideal for cutting 1/16-inch tubing.
- Cuts quickly, reducing distortion.
- Cuts clean, eliminating need for deburring.



Description

1/16" Tubing Cutting Pliers

qty.

ea.

cat.#

20193

price

Clean-Cut Tubing Cutter

- Burr-free, perpendicular cuts that will not distort the tubing OD or close the ID.
- Use on PEEK™, Teflon®, Tefzel®, other polymeric tubing.

25069

Description

Clean-Cut Tubing Cutter

Replacement Blade for Clean-Cut Cutter

qty.

ea.

ea.

cat.#

25069

25070

price



*Not for use with stainless steel tubing.



SSI TC-20 Tube Cutting Machine

- Cuts $\frac{1}{16}$ ", $\frac{1}{8}$ ", or $\frac{1}{4}$ " tubing with inside diameter as small as 0.008".
- Electrically operated bench-top model.
- Handy dressing tool on the swing arm removes burrs and reams tubing.
- Voltage selectable 110–120/220–240 volts, 50–60 Hz.*

| Description | qty. | cat.# | price |
|--------------------------------|-------|-------|-------|
| SSI Tubing Cutter Machine | ea. | 23029 | |
| SSI Replacement Cutting Wheels | 3-pk. | 23030 | |

*Unit shipped set for 110–120 operating voltage. Switch to 220–240 volts by using alternate fuse and power cord (included).

Tubing Burring & Reaming Tool

Removes burrs and reams tubing.



20188

| Description | Size | qty. | cat.# | price |
|----------------------|--------------------------------|------|-------|-------|
| Tubing Dressing Tool | 1/16" | ea. | 20188 | |
| Replacement Insert | for 1/16" Tubing Dressing Tool | ea. | 20189 | |
| Tubing Dressing Tool | 1/8" | ea. | 20190 | |
| Replacement Insert | for 1/8" Tubing Dressing Tool | ea. | 20191 | |

Tubing Reamer

- Removes burrs from stainless steel tubing.
- For 1/4- or 1/8-inch tubing.
- Nonslip safety design.



Description

Tubing Reamer

| qty. | cat.# | price |
|------|-------|-------|
| ea. | 20134 | |

Rethreading Tool

- Repair worn or damaged threads.
- Multiple uses (fittings, injectors, detectors, etc.)
- Built-in guide to prevent cross-threading.



Description

Rethreading Tool for 1/16" compression fitting

Rethreading Tool for 1/8" compression fitting

Rethreading Tool for 1/4" compression fittings
(Agilent split/splitless injection ports)

Rethreading Tool

for 7/16" compression fitting (Varian injection ports)

Rethreading Tool

for 1/4" Varian-style capillary column fittings

| qty. | cat.# | price |
|------|-------|-------|
| ea. | 23016 | |
| ea. | 23017 | |
| ea. | 23018 | |
| ea. | 23019 | |
| ea. | 21893 | |

Pocket Magnifier

- Small and easy to handle.
- 10x magnification makes it easy to see the column end to verify a square cut.



Description

Pocket Magnifier

| qty. | cat.# | price |
|------|-------|-------|
| ea. | 20124 | |

Flashlight with Swivel Head

- Ideal for tight spaces—like inside a GC oven.
- Uses two AA batteries (included).



Description

Flashlight with Swivel Head

| qty. | cat.# | price |
|------|-------|-------|
| ea. | 22187 | |

Flexible Inspection Light

- Inspect inside surfaces of sample cylinders or other chambers.
- 14" reach.
- 100,000-hour LED life.



Description

Flexible Inspection Light

| qty. | cat.# | price |
|------|-------|-------|
| ea. | 22627 | |