


**Figure 2** OPN on Res-Sil™ C packing has unique selectivity for *cis*-2-butene and 1,3-butadiene.




**Darrel Zezzo**  
Ohio Valley Sales  
Representative  
1+ year of service!

also available

Res-Sil™ packing materials.  
See page 121.

**Figure 3** Refinery gas calibration standard on a Restek refinery gas packed column set.



2abc Refinery Gas Column Set (cat.# 88000-875) (3 column set)  
 Oven temp.: 60°C  
 Inj. temp.: 150°C  
 Det. temp.: 150°C FID  
 Flow rate: 30mL/min., helium  
 Sample: refinery gas  
 Sample size: 1cc