

## Custom Packed Columns

### To order, specify the following:

- 1) column dimensions (length, ID) and tubing material
- 2) packing description (percent coating and phase, support mesh size, and treatment)
- 3) column configuration (instrument manufacturer, model number, on-column injection or not) and with or without nuts and ferrules

**Ordering Example:** (6' x 1/8") (stainless steel) (3%) (Rtx®-1) (Silcoport™ P) (80/100) (Agilent 6890) (on-column injection) (fittings kit).

Please use the custom order form on page 128.

### ordering note

For international pricing on custom packed or micropacked columns, please contact your Restek representative.

## Custom Micropacked Columns

### To order, contact your Restek representative and specify the following:

- 1) physical dimensions (length, OD, ID, and tubing material)
- 2) packing description (percent coating and phase, support mesh size, and treatment)
- 3) installation kit (see page 119)

**Ordering Example:** (2m x 1/16" OD x 1.00mm ID) (Siltek®-treated tubing) (5%) (Carbowax® 20M) (CarboBlack B) (80/120) (installation kit for valve applications, cat. #21065)

Please use the custom order form on page 128.