

## Capillary Column Nuts for Agilent 5890/6890/7890 GCs

Available in brass or stainless steel.



Description	Similar to Agilent part #	qty.	cat.#
<b>For use with "compact" Agilent-style ferrules.</b>			
Stainless Steel Capillary Column Nut	5181-8830	2-pk.	21884
Brass Capillary Column Nut	5181-8830	2-pk.	21878
<b>For use with standard 1/16"-type ferrules.*</b>			
Stainless Steel Capillary Column Nut	05921-21170	2-pk.	20883
Brass Capillary Column Nut	05921-21170	2-pk.	21879



for use with "compact" Agilent-style ferrules



for use with standard ferrules

\*Designed to fit a wider variety of 1/16" ferrules

## Capillary Column Nuts for Agilent 6850 GCs

- Meet or exceed original equipment performance.
- High-quality stainless steel.

Description	Similar to Agilent part #	qty.	cat.#
Capillary Column Nuts for Agilent 6850 GCs	5183-4732	2-pk.	21882



## Finger-Tight Capillary Column Nuts

- Rapidly tighten without wrenches; avoid stripped threads.
- Either version can be used with 0.25, 0.32, or 0.53mm ID columns.

Description	qty.	cat.#
<b>For use with "compact" Agilent-style ferrules.</b>		
Finger-Tight Capillary Column Nut	ea.	21311
<b>For use with standard ferrules.</b>		
Finger-Tight Capillary Column Nut	ea.	21312



cat.# 21312

## Finger-Tight Capillary Column Nuts


- Allow wrench-free column installations.
- Work with standard or compact (Agilent-style) ferrules.
- Made from high-quality stainless steel.



Description	Similar to Agilent part #	qty.	cat.#
<b>For use with "compact" Agilent-style ferrules.</b>			
A) Finger-Tight Capillary Column Nuts	5020-8293	2-pk.	21040
	5020-8292		
<b>For use with standard 1/16"-type ferrules.</b>			
B) Finger-Tight Capillary Column Nuts	5020-8293	2-pk.	21041
	5020-8292		
<b>For use with standard 1/16" compression fittings</b>			
C) Finger-Tight Capillary Column Nuts		2-pk.	21042

also available

**Ferrules**  
See pages 217-219 for our complete listing of ferrules.



## Capillary Installation Gauge

- Seats graphite\* ferrules onto column for consistent installations.
- Prevents crushed column ends.
- Made from high-quality stainless steel.



Thread...



...tighten...



...done.



A  
B

Easily seat ferrules for consistent installations!

Description	qty.	cat.#
A) Capillary Installation Gauge for Agilent-style fittings (compact ferrules)	ea.	21034
B) Capillary Installation Gauge for 1/16" fittings (1/16" ferrules)	ea.	21399

\*Recommended for use with graphite ferrules