


SureFit Split/Splitless Injection Ports for Agilent 5890 GCs

restek
exclusive!

GC ACCESSORIES | INSTRUMENT SUPPLIES

Uses Restek Enhanced Base Fitting (cat.# 22323)

Uses Manufacturer's-Style Inlet Seal and Reducing Nut (cat.# 22326)



Donna Lidgett
GC Accessories Product
Marketing Manager
21+ years of service!

SureFit Split/Splitless Injection Ports for Agilent 5890 GCs

Uses Restek Enhanced Base Fitting

Description	qty.	cat.#
E) Injection Port Assembly Kit includes: * Siltek® split/splitless weldment, Siltek® shell weldment, base fitting, Siltek® base screw, septum nut, 1/16" and 1/4" stainless steel nuts, 1/4" graphite ferrule	kit	22323
A) Siltek® Split/Splitless Weldment for Restek 5890 Injection Port	ea.	22324
B) Siltek® Shell Weldment for Restek 5890 Injection Port	ea.	22325
Ferrules for split liners (6.3mm OD): 1/4" dual-taper VespeL®/graphite	10-pk.	20290
Ferrules for splitless liners (6.5mm OD): 1/4" dual-taper VespeL®/graphite	10-pk.	20291
C) Siltek® Base Screws for Restek Injection Port for 5890/6890	10-pk. 50-pk.	21631 21632
Gold-Plated Base Screws for Restek Injection Port for 5890/6890	2-pk. 10-pk.	21629 21630
Autosampler & PTV Septum Nut (for 23-gauge needles)	ea.	20631
D) Base Fitting for Restek Injection Port for 5890/6890	ea.	21626

Uses Manufacturer's-Style Inlet Seal and Reducing Nut

Description	qty.	cat.#
F) Injection Port Assembly Kit includes: * Siltek® split/splitless weldment, Siltek® shell weldment, septum nut, Siltek® Dual VespeL® Ring Inlet Seal, Agilent-style stainless steel reducing nut, Agilent-style stainless steel capillary nut	kit	22326
A) Siltek® Split/Splitless Weldment for Restek 5890 Injection Port	ea.	22324
B) Siltek® Split/Splitless Shell Weldment for Restek Injection Port for 5890/6890 (uses standard Agilent inlet seals)	ea.	22327
Ferrules for split liners (6.3mm OD): 1/4" dual-taper VespeL®/graphite	10-pk.	20290
Ferrules for splitless liners (6.5mm OD): 1/4" dual-taper VespeL®/graphite	10-pk.	20291
Autosampler & PTV Septum Nut (for 23-gauge needles)	ea.	20631
D) Reducing Nut	ea.	22078
E) Stainless Steel Capillary Column Nut	2-pk.	21884
C) Siltek® 0.8mm ID Dual VespeL® Ring Inlet Seal**	2-pk. 10-pk.	21242 21243



please note

Special dual-taper, 1/4-inch ID VespeL®/graphite ferrules fit over split liners (6.3mm OD) (cat.# 20290). Ferrules with slightly enlarged ID, required for splitless liners (6.5mm OD) also available (cat.# 20291).

*Injection port assembly kit does not include inlet liner, 1/16" capillary ferrule, or split/splitless liner ferrules. Order separately.

**For more inlet seals, refer to pages 150-151.