

**Oven Flapper Assemblies**

- Ensure even and consistent column heating, run after run, by replacing the oven flapper assembly when the insulation becomes damaged or worn.
- Replacement Insulation Gasket and Hardware Kit includes the insulation gasket, 2 screws and 2 nuts.
- Meets or exceeds manufacturer's performance.
- Insulation is replaceable, using Restek's Oven Flapper Assemblies.

Description	Similar to Agilent part #	qty.	cat.#
Oven Flapper Assembly for Agilent 5890 GCs	05890-80560	ea.	22994
Oven Flapper Assembly for Agilent 6890 GCs	G1530-80560	ea.	22995
Replacement Insulation Gasket and Hardware Kit*		kit	22996

\*Fits Restek's Oven Flapper Assembly only.

**Air Diverter for Agilent GCs**

- Divert GC exhaust heat away from the lab bench.
- Meets or exceeds manufacturer's performance.
- Easy to install—no tools required.

Description	Similar to Agilent part #	qty.	cat.#
Air Diverter for Agilent 5890/6890 GCs	19247-60510 G1530-80650	ea.	22076

**Make Life Easier (MLE) Capillary Tool Kit for Agilent GCs****Includes:**

- capillary installation gauge for Agilent GCs
- injector wrench for Agilent GCs
- septum nut removal tool
- 1/8", 3/16", and 1/4" nylon brushes
- 1/4", 3/8", and 3/16" stainless steel wire tube brushes
- stainless steel surface brush
- 6 stainless steel jet reamers (0.25–0.65mm OD)
- 3/4" x 5/16" open end wrench
- 3/8" x 7/16" open end wrench
- 7/16" x 1/2" open end wrench
- 1/2" x 9/16" open end wrench
- rubber-tipped slide-lock tweezers
- scoring wafers with handles
- inlet liner removal tool
- septum puller
- mini wool puller/insertor tool
- 4-inch tapered needle file
- swivel head flashlight
- mini hand drill set
- 15cm compact steel ruler
- pocket magnifier
- high temperature string (1 meter)
- pipe cleaner (12-inch)
- cotton tip swabs (pk. of 25)



Description	qty.	cat.#
MLE Capillary Tool Kit for Agilent GCs	kit	22186

Everything you need  
in one complete kit!