

## Direct Injection (DI) Inserts and Liner Adaptor for 1/4-Inch Packed Column Injection Ports (for 0.25, 0.32, & 0.53mm ID columns)

### DI Glass Inserts for Agilent 5890 Packed Column GC

- Tolerances closely controlled.
- Can be removed from the septum nut weldment.



Description	Similar to Agilent part #	qty.	cat.#
DI Glass Inserts for Agilent 5890 Packed Column GC (1.7mm ID x 3.0mm OD x 93mm)	5181-3382 5080-8732	5-pk.	20967

All liners are  
**100%**  
deactivated

All liners are shipped intermediate polarity (IP) deactivated unless otherwise requested.

**SILTEK**<sup>®</sup>

version available.

### DI Uniliner<sup>®</sup> Liners for Agilent 5890 Packed Column GC

- Press-Tight<sup>®</sup> taper forms dead-volume-free connection to column.
- Minimizes solvent and peak tailing.
- Use with 0.25, 0.32, and 0.53mm ID capillary columns.
- Can be removed from the septum nut weldment.



Description	Fits same liner adaptor as Agilent part #	qty.	cat.#
DI Uniliner <sup>®</sup> Liner for Agilent 5890 Packed Column GC (1.7mm ID x 3.0mm OD x 93mm)	5181-3382 5080-8732	ea.	20964
DI Uniliner <sup>®</sup> Liners for Agilent 5890 Packed Column GC (1.7mm ID x 3.0mm OD x 93mm)	5181-3382 5080-8732	5-pk.	20965
DI Uniliner <sup>®</sup> Liners for Agilent 5890 Packed Column GC (1.7mm ID x 3.0mm OD x 93mm)	5181-3382 5080-8732	25-pk.	20966

### DI Liner Adaptor for Agilent 5890 Packed Column GC

- Uses standard 1/16-inch capillary nut and ferrules.
- Convenient wrench pad at base.
- Includes 1/4-inch graphite ferrule and stainless steel nut.
- Use with Agilent or Restek DI glass inserts or Restek DI Uniliner<sup>®</sup> liners for an Agilent 5890 packed column GC.

Description	Similar to Agilent part #	qty.	cat.#
DI Liner Adaptor for Agilent 5890 Packed Column GC	19244-80540	ea.	21303



also **available**

#### Looking for autosampler vials?

Restek carries a full line of autosampler vials for all of Agilent's autosamplers. See [pages 259 to 275](#) or visit [www.restek.com](http://www.restek.com) and search on autosampler vial.