

Replacement Cables for Agilent GCs, Integrators, & Autosamplers

- Priced less than original equipment.
- Instructions and wiring diagrams included.
- Manufactured with only the highest-quality components.

did you **know?**

Restek provides pin schematics.

Connect an Agilent 5890 GC to an Agilent integrator (for second Inet integrator).

| Description | Similar to Agilent part # | length | qty. | cat.# |
|--------------------|---------------------------|--------|-------|-------|
| Replacement Cable | 35900-60610 | 6 ft. | ea. | 20650 |
| Replacement Cables | 35900-60610 | 6 ft. | 2-pk. | 20651 |



Connect an Agilent 5890 GC to an Agilent 3396 integrator to enable remote starts (non-Inet connection from GC to integrator).

| Description | Similar to Agilent part # | length | qty. | cat.# |
|-------------------|---------------------------|--------|------|-------|
| Replacement Cable | 03394-60560 | 6 ft. | ea. | 20654 |



Connect an Agilent 3396 integrator to remote-start another piece of equipment or to start the Agilent 3396 integrator from that other piece of equipment.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|--------------------|---------------------------|--------|-------|-------|
| Replacement Cable | 03394-60540 | 6 ft. | ea. | 20655 |
| Replacement Cables | 03394-60540 | 6 ft. | 2-pk. | 20656 |



Connect an Agilent 5890 GC to a non-Agilent integrator or standard strip chart recorder.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|--------------------|---------------------------|--------|-------|-------|
| Replacement Cable | 05890-60800 | 6 ft. | ea. | 20652 |
| Replacement Cables | 05890-60800 | 6 ft. | 2-pk. | 20653 |



Connect an Agilent 5890 GC to remote-start another piece of equipment or to start the Agilent 5890 GC from that piece of equipment.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|-------------------|---------------------------|--------|------|-------|
| Replacement Cable | 05890-61080 | 6 ft. | ea. | 20657 |



Connect an Agilent 3396 integrator to a non-Agilent GC.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|--------------------|---------------------------|--------|-------|-------|
| Replacement Cable | 35900-60630 | 6 ft. | ea. | 20658 |
| Replacement Cables | 35900-60630 | 6 ft. | 2-pk. | 20659 |



Connect an Agilent 5890 GC to an Agilent 35900C Interface or to an Agilent 7673A Autosampler to enable remote starts.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|--------------------|---------------------------|--------|-------|-------|
| Replacement Cable | 35900-60700 | 6 ft. | ea. | 20660 |
| Replacement Cables | 35900-60700 | 6 ft. | 2-pk. | 20661 |



Automatic liquid sampler stop/start cable for Agilent GCs.

| Description | Similar to Agilent part # | length | qty. | cat.# |
|-------------------|---------------------------|--------|------|-------|
| Replacement Cable | G1530-60930 | 7 ft. | ea. | 20929 |



Connect an Agilent 6890 GC to a PC 9F/9F (RS232).

| Description | Similar to Agilent part # | length | qty. | cat.# |
|-------------------|---------------------------|--------|------|-------|
| Replacement Cable | G1530-60600 | 6 ft. | ea. | 22694 |

