

**FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL**

- Made of high-quality stainless steel.
- Meet or exceed manufacturer's performance.

Description	Similar to PE part #	qty.	cat.#
<b>A) For use with PE style capillary nuts</b>			
FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL	N6120020	ea.	22608
<b>B) For use with 1/16" compression style nuts</b>			
FID Capillary Column Adaptor for PerkinElmer Auto SYS™ XL	—	ea.	22609

**FID Replacement Parts for PerkinElmer Auto SYS™ XL and Clarus 500**

- Made of high-quality stainless steel.
- Meet or exceed manufacturer's performance.



Description	Similar to PE part #	qty.	cat.#
FID Jet for PerkinElmer Auto SYS™ XL and Clarus 500	N6100361	ea.	23038
Auto-Ignite FID Replacement Part Kit for PerkinElmer Auto SYS™ XL and Clarus 500	N6103167	kit	23061
	N6103175		
	N6101085		
	N6001204		
09912223			
Nozzle Insulator for PerkinElmer Auto SYS™ XL and Clarus 500	N6101085	ea.	23062
FID Body for PerkinElmer Auto SYS™ XL and Clarus 500	N6100364	ea.	23063

**Make Life Easier (MLE) Capillary Tool Kit for PerkinElmer GCs****Includes:**

- 1/8", 3/16", and 1/4" nylon brushes
- 1/4", 3/8", and 3/16" stainless steel wire tube brushes
- stainless steel surface brush
- 6 stainless steel jet reamers (0.25–0.65mm OD)
- 1/4" x 5/16" open end wrench
- 3/8" x 7/16" open end wrench
- 7/16" x 1/2" open end wrench
- 1/2" x 9/16" open end wrench
- rubber-tipped slide-lock tweezers
- scoring wafers with handles
- inlet liner removal tool
- septum puller
- mini wool puller/insertor tool
- 4-inch tapered needle file
- swivel head flashlight
- mini hand drill set
- 15cm compact steel ruler
- pocket magnifier
- high temperature string (1 meter)
- pipe cleaner (12-inch)
- cotton tip swabs (pk. of 25)

Everything you need  
in one complete kit!



Description	qty.	cat.#
MLE Capillary Tool Kit for PerkinElmer GCs	kit	22185