

## Septa for Thermo Scientific GCs

- Preconditioned and ready to use.
- Packaged in precleaned glass jars.

Now precision  
molded!

Septum Diameter	25-pk.	50-pk.	100-pk.
<b>Thermolite® Septa (usable to 340°C inlet temp.)</b>			
9mm	27132	27133	27134
9.5mm (9/16")	27135	27136	27137
17mm	27150	27151	27152



Thermolite® Septa

<b>IceBlue™ Septa (usable to 250°C inlet temp.)</b>			
9mm	—	27156	27157
9.5mm (9/16")	—	27158	27159
17mm	—	27168	27169



IceBlue® Septa

<b>BTO® Septa (usable to 400°C inlet temp.)</b>			
9mm CenterGuide™	—	27104	27105
9.5mm (9/16")	—	27106	27107
17mm CenterGuide™	—	27116	27117



BTO® Septa

## HANDY septum size chart

Instrument	Septum Diameter (mm)
<b>Thermo Scientific</b>	
TRACE™ GC, GCQ w/TRACE™, PTV, 8000 series	17
<b>Finnigan (TMQ)</b>	
GC 9001, GCQ, QCC™, TRACE™ 2000	9.5

## Septa Holder Kits for Thermo Scientific TRACE™ and Focus SSL

- Includes septum support and holder.
- Made from high-quality stainless steel.

Description	Similar to TS part #	qty.	cat.#
A) Septa Holder for Thermo Scientific TRACE™ and Focus SSL	233 030 15		
B) Silcosteel®-AC Treated Septa Holder for Thermo Scientific TRACE™ and Focus SSL	350 054 335	kit	21299
	233 030 15		
	350 054 33	kit	24972

**new!**



## Septum Cap for Split/Splitless Injector on Thermo Scientific TRACE™ and Focus SSL

Description	Similar to TS part #	qty.	cat.#
Septum Cap for Split/Splitless Injector on Thermo Scientific TRACE™ and Focus SSL	350 01 050	ea.	24971

**new!**



## Silver Seals for Thermo Scientific TRACE™ and Focus SSL

Description	Similar to TS part #	qty.	cat.#
A) Silver Seal for Split/Splitless Injector	29033629	10-pk.	23057
A) Silver Seal for Split/Splitless Injector	29033629	20-pk.	23058
B) Silver Seal for Capillary Detector Base Body	29037100	10-pk.	23059
B) Silver Seal for Capillary Detector Base Body	29037100	20-pk.	23060

