

## ELCD/PFPD Replacement Parts



**Don McCandless**  
Head Coach  
18+ years of service!

### ELCD Nickel Reaction Tubes

- Pretreated for maximum sensitivity.
- Quality-controlled for reliability.
- Available for many popular models.



ELCD Model #	To replace these instrument part numbers:					Order these Restek part numbers:	
	Tremetrics	Varian	PerkinElmer	Shimadzu	O.I. Analytical	qty.	cat.#
Hall 1000	117459-0003	00-997625-12	N660-1072	220-90435-00	—	2-pk.	21581
O.I. 4420	—	—	—	—	260323	2-pk.	21582

### Replacement Accessories for Hall 1000

#### ELCD Nickel Reaction Tube Nut

High-quality stainless steel ELCD nut mounts nickel reaction tube into ELCD.



Description	qty.	cat.#
ELCD Nickel Reaction Tube Nut	2-pk.	21584

#### 1/16-Inch Vespel®/Graphite Sealing Ring

Installs onto the nickel reaction tube after the screw. Easily compresses on the reaction tube to provide a leak-tight seal and prevent detector oxidation.



Description	qty.	cat.#
1/16-Inch Vespel®/Graphite Sealing Ring	2-pk.	21583

#### Cleaned Teflon® Transfer Lines for ELCDs

We stringently clean our ELCD Teflon® transfer lines with an HCl solution to remove any contaminants, then rinse with methanol. Convenient 6.5-inch pre-cut pieces directly interface the nickel reaction tube and conductivity cell in Tracor, Tremetrics, O.I., and many other ELCDs.



Description	qty.	cat.#
Teflon® Transfer Lines for ELCDs, 1/16" OD x 0.020" ID (five 6.5-inch lines)	5-pk.	20121