

Make Life Easier (MLE) Capillary Tool Kits

Everything you need in one complete kit!

All kits include these components:

- $\frac{1}{8}$ " $\frac{3}{16}$ " $\frac{1}{4}$ " nylon brushes
- $\frac{1}{4}$ " $\frac{3}{8}$ " $\frac{3}{16}$ " stainless steel wire tube brushes
- stainless steel surface brush
- 6 stainless steel jet reamers (0.25–0.65mm OD)
- $\frac{3}{4}$ " x $\frac{7}{16}$ " open end wrench
- $\frac{3}{8}$ " x $\frac{7}{16}$ " open end wrench
- rubber-tipped slide-lock tweezers
- scoring wafers with handles
- inlet liner removal tool
- septum puller
- mini wool puller/insertor tool
- 4-inch tapered needle file
- swivel head flashlight
- mini hand drill set
- 15cm compact steel ruler
- pocket magnifier
- high temperature string (1 meter)
- pipe cleaner (12-inch)
- cotton tip swabs (pk. of 25)

MLE Capillary Tool Kit for Agilent GCs (cat.# 22186) also includes:

- capillary installation gauge for Agilent GCs
- injector wrench for Agilent GCs
- septum nut removal tool
- $\frac{7}{16}$ " x $\frac{1}{2}$ " open end wrench
- $\frac{1}{2}$ " x $\frac{7}{16}$ " open end wrench

MLE Capillary Tool Kit for PerkinElmer GCs (cat.# 22185) also includes:

- $\frac{7}{16}$ " x $\frac{1}{2}$ " open end wrench
- $\frac{1}{2}$ " x $\frac{7}{16}$ " open end wrench

MLE Capillary Tool Kit for Shimadzu GCs (cat.# 22182) also includes:

- capillary installation gauge for Shimadzu GCs
- injector wrench for Shimadzu GCs
- 6mm x 7mm open end wrench
- 8mm x 10mm open end wrench
- 16mm x 17mm open end wrench

MLE Capillary Tool Kit for Thermo Scientific GCs (cat.# 22183) also includes:

- capillary installation gauge for Thermo Scientific GCs
- liner cap removing tool for Thermo Scientific GCs
- 6mm x 7mm open end wrench
- 8mm x 10mm open end wrench
- 16mm x 17mm open end wrench

MLE Capillary Tool Kit for Varian GCs (cat.# 22184) also includes:

- capillary installation gauge for Varian GCs
- $\frac{7}{16}$ " x $\frac{1}{2}$ " open end wrench
- $\frac{1}{2}$ " x $\frac{7}{16}$ " open end wrench



For Agilent GCs
(cat.# 22186)



For PerkinElmer GCs
(cat.# 22185)



For Shimadzu GCs
(cat.# 22182)



For Thermo Scientific GCs
(cat.# 22183)



For Varian GCs
(cat.# 22184)

Description	qty.	cat.#
MLE Capillary Tool Kit for Agilent GCs	kit	22186
MLE Capillary Tool Kit for PerkinElmer GCs	kit	22185
MLE Capillary Tool Kit for Shimadzu GCs	kit	22182
MLE Capillary Tool Kit for Thermo Scientific GCs	kit	22183
MLE Capillary Tool Kit for Varian GCs	kit	22184

did you know?

Make Life Easier!

MLE Tool Kits conveniently provide the tools that make it easier to install and maintain capillary columns!