


Super-Clean™ Gas Filters and Base Plates



All traps measure:
10⁵/₈" x 1³/₄" (27 x 4.4 cm)
Each base plate unit
measures:
4" x 4" x 1⁷/₈"
(10.2 x 10.2 x 4.8 cm)

Super-Clean™ Gas Filter and Base Plate Kits

Description	qty.	cat.#
Carrier Gas Cleaning Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Triple Gas Filter)	kit	22019
Fuel Gas Purification Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and hydrocarbon/moisture Fuel Gas Filter)	kit	22021

Replacement Gas Filters

Description	qty.	cat.#
Replacement Triple Gas Filter (removes oxygen, moisture and hydrocarbons)	ea.	22020
Replacement Fuel Gas Filter (removes moisture and hydrocarbons)	ea.	22022

Gas Filter Bundle Kit

- Ideal for use in combination with 3-position base plate—purchase separately.

Description	qty.	cat.#
Gas Filter Bundle Kit	kit	22031

Super-Clean™ Ultra-High Capacity Gas Filters

Description	qty.	cat.#
Ultra-High Capacity Hydrocarbon Filter	ea.	22030
Ultra-High Capacity Moisture Filter	ea.	22028
Ultra-High Capacity Oxygen Filter	ea.	22029

Helium-Specific Super-Clean™ Gas Filter and Kit

- Specifically designed for purification of helium in GC/MS systems.
- Traps are packed and conditioned using helium.
- Uses standard single-position base plate.

Description	qty.	cat.#
Helium-Specific Carrier Gas Cleaning Kit (includes mounting base plate, 1/8" inlet/outlet fittings, and oxygen/moisture/hydrocarbon Helium-Specific Filter)	kit	21983
Replacement Helium-Specific Gas Filter (removes oxygen, moisture and hydrocarbons)	ea.	21982

Filter Base Plates

- Standard base plate fittings are 1/8". To adapt to 1/4", order 1/8" to 1/4" tube-end unions.
- Base plates fit all Super-Clean™ gas filters listed above.

Description	qty.	new!	
		Brass cat.#	Stainless Steel cat.#
Single-Position Filter Base Plate	ea.	22025	22344
2-Position Filter Base Plate	ea.	22026	22345
3-Position Filter Base Plate	ea.	22027	22346

Wall Mounting Bracket

- Base plates can be mounted by using screws and the mounting holes on the base plate, or by using this optional wall mounting bracket.

Description	qty.	cat.#
Wall Mounting Bracket for Super-Clean™ Base Plates	ea.	21984

Replacement O-Rings for Cartridge Base Plates

- Pack includes 10 large O-rings and 10 small O-rings.

Description	qty.	cat.#
Replacement O-Rings for Cartridge Base Plates	20-pk.	22023

1/8-Inch to 1/4-Inch Tube-End Unions

- To adapt 1/8" Super-Clean™ base plate fittings to 1/4", use 1/8" to 1/4" tube-end unions.

Description	qty.	Stainless Steel	
		Brass cat.#	cat.#
1/8" to 1/4" Tube-End Unions	5-pk.	21833	21933