

## CGA Fittings

CGA-specified nuts and nipples with internal frit, 1/4-inch NPT nickel-plated brass.

Description	qty.	cat.#
CGA 580 Fitting, (N <sub>2</sub> , He, Ar)	ea.	21336
CGA 350 Fitting, (H <sub>2</sub> , P <sub>2</sub> )	ea.	21337
CGA 590 Fitting, (Air)	ea.	21338



## Flexible Stainless Steel Hoses

Description	qty.	cat.#
Flexible Stainless Steel Hose, 36", 1/4" Female NPT Fittings	ea.	21339
Flexible Stainless Steel Hose, 18", 1/4" Female NPT Fittings	ea.	21340
Flexible SS Hose, 36", with Stainless Steel CGA 580	ea.	21344
Flexible SS Hose, 36", with Stainless Steel CGA 350	ea.	21345
Flexible SS Hose, 36", with Stainless Steel CGA 590	ea.	21346



## Flammable Gas Flash Arrestor—Factory Mutual Approved\*

- Gas flow shuts off in the event of a flashback.
- Flame extinguished—flame front prevented from reaching the gas supply.
- No gas flow restriction under normal operating conditions.

Description	qty.	cat.#
Flammable Gas Flash Arrestor, Brass Body	ea.	21334

\*Approved for brass body servicing hydrogen, acetylene, propane, or natural gas only.



## Backpressure Gas Regulator

Capillary GC inlet systems have backpressure regulators to maintain a constant upstream pressure and rapidly respond to catastrophic leaks. The 0–60psig (0–414kPa) operating range is sufficient to operate a 105m, 0.25mm ID column at its optimum flow rate.

Description	qty.	cat.#
Backpressure Gas Regulator	ea.	20635



## MINICYL Regulator

This compact general purpose regulator has many laboratory applications including air-drying glassware, sparging or evaporating solutions, and controlling pneumatic valves. It is constructed of lightweight aluminum with an elastomer diaphragm. Includes a 0–60psig (0–414kPa) gauge and either 1/8- or 1/4-inch tube fittings.

Description	qty.	cat.#
MINICYL Regulator, 1/8" Fittings	ea.	20610
MINICYL Regulator, 1/4" Fittings	ea.	20611



## Cylinder Valve Wrench

This specially-designed wrench enables easy opening of cylinder valves that are fitted with a hand wheel. It is also suitable for removing difficult cylinder caps.

Description	qty.	cat.#
Cylinder Valve Wrench	ea.	21321



## Universal Cylinder Wrench

Use this versatile wrench for tightening gauges and gas regulator CGA fittings to cylinder outlets and pipe thread connections.

Description	qty.	cat.#
Universal Cylinder Wrench	ea.	21322



## Wall-Mount Gas Cylinder Holder

Prevent serious injuries! This holder is designed to prevent free-standing gas cylinders from tipping over and injuring personnel. The cast aluminum holder can be secured to a wall or the side of a work bench. The mount will secure a cylinder from 4–15 inches in diameter.

Description	qty.	cat.#
Wall-Mount Gas Cylinder Holder	ea.	21333

