

Custom Phases

Over the years, we have developed numerous proprietary phases for specific separations. We have more than 35 custom stationary phases in stock. Also note that some of our “common” phases are available in variations for specific purposes. If you don't see it in this catalog, call us.



Paul Silvis
 Founder, Inventor, Coach
 23+ years of service!

Custom Phases Currently in Stock*

Apiezon® L	DEGS (G4)	OV®-351 (G35)	Squalane
Butane 1,4 diol succinate	DEXSIL 300	Rt™-2330 (G8)	Squalene
Carbowax® 300	DIIDP (G24)	Rt™-2340 (G5)	Superox® 20M
Carbowax® 600	FFAP	SE®-30 (G1,G2)	TCEP
Carbowax®1000 (G14)**	OS®-124	SE®-54 (G36)	Tricresyl phosphate
Carbowax® 8000 (G47)	OS®-138	Silar® 5CP	UCON® 50HB
Carbowax® 20M (G16)	OV®-202	Silar® 9CP (G8)	XE®-60
DC-550 (G28)	OV®-275	Silar® 10CP	XF®-1150

*Other phases available on request.

**Codes in parentheses are USP designations—see page 125.

Custom Column Dimensions

Need a specific film thickness? Locked into a method that uses a 42-meter column? If you need it, we can make it. We have shipped lengths from 3 to 150 meters, and film thicknesses from 0.1µm to 10.0µm.

Inverse Gas Chromatography Columns

We have extensive experience in coating capillary columns with a wide variety of polymers. If you need an IGC column, contact Technical Service at 800-356-1688 or 814-353-1300, ext. 4, or contact your Restek representative. We'll assist you in coating a column with your polymer.

PIE service

Plus 1

Exceeding your expectations in everything we do.

Innovation

Turning Visions into Reality®.

Execution

On-time delivery of products and services.

Restek's vision is to be the company that chromatographers trust by providing the highest quality, most innovative products and services throughout the world.

We will soon reach our goal of 100% employee ownership. As owners, our success depends on your success.

