

Convenience Kits: Vials, Caps, & Septa

Vials packaged in a clear-lid tray. Caps with septa packaged in a plastic bag.

Description	100-pk.	1000-pk.
2.0mL Clear Vial, deactivated , PTFE/Natural Rubber Silver Seal	24671	24672
2.0mL Amber Vial, deactivated , PTFE/Natural Rubber Silver Seal	24673	24674
2.0mL Clear Vial, untreated, PTFE/Natural Rubber Silver Seal	21196	21197
2.0mL Amber Vial, untreated, PTFE/Natural Rubber Silver Seal	21198	21199
2.0mL Clear Vial, untreated, PTFE/Silicone Silver Seal	24646	24647
2.0mL Amber Vial, untreated, PTFE/Silicone Silver Seal	24648	24649

**Limited Volume Inserts for 2.0mL Crimp-Top & Short-Cap, Screw-Thread Vials**

Description	100-pk.	1000-pk.
50µL Glass, Big Mouth w/ Bottom Spring	24513	21782
250µL Glass, Big Mouth Insert w/ Bottom Spring	21776	21777
250µL Glass, Big Mouth Insert w/ Glass Flange (Step™ Design)*	24516	21779
350µL Glass, Flat Bottom Insert	21780	24517
350µL Glass, Flat Bottom Insert w/ ID Ring	24692	24693
250µL Polypropylene, Big Mouth, Bottom Spring	24518	—
250µL Polypropylene, Big Mouth, Top Flange	24519	—
250µL Polypropylene, Big Mouth, No Spring	24520	—

*Not to be used with 9mm screw-thread vials.

**Limited Volume Vials, 11mm Crimp-Top**

Description	100-pk.	1000-pk.
100µL Polypropylene	24651	24652
100µL Glass Insert/Clear Plastic Vial	24653	—

**Versa Vial™ Vials, Glass, 2.0mL, 12 x 32mm**

- Increased neck length for accommodating fluctuations in robotic arm settings.
- Eliminate cracked necks, breakage and dangerous sharp edges, using soft Versa Vial™ PTFE/silicone plugs.
- Largest target area available, ideal for sampling and pipet use.
- Vials available with or without graduated marking spots.
- Patent pending.

Description	100-pk.	1000-pk.
2.0mL Clear Glass Vial	22701	22702
2.0mL Clear Glass Vial w/Graduated Marking Spot	22703	22704
2.0mL Amber Glass Vial	22697	22698
2.0mL Amber Glass Vial w/Graduated Marking Spot	22699	22700

**Versa Vial™ Vials, Polypropylene, 1.5mL, 12 x 32mm**

Description	100-pk.	1000-pk.
1.5mL Polypropylene Vial	22705	22706

**Versa Vial™ Inserts**

Description	100-pk.	1000-pk.
250µL Glass Insert w/Glass Flange	22707	22708
250µL Polypropylene w/Flange	22709	22710
250µL Glass Flat Bottom Insert	22711	22712

**Versa Vial™ Plugs**

Description	100-pk.	1000-pk.
12mm, Green Polyethylene Plug	22713	22714
12mm, Grey Chlorobutyl/Siliconized Plug	22715	22716
12mm, White PTFE/Silicone Plug	22717	22718
12mm, White PTFE/Silicone Plug w/Slit*	22719	22720



*Eliminates needle damage to plug.

2.0mL, 11mm Poly Crimp™ Vials



2.0mL, 11mm Poly Crimp™ Seal: Snap-on or Crimp

The patented Poly Crimp™ Seal is versatile, working as either a snap-on or crimp-top cap. Simply use a standard crimping tool to secure the Poly Crimp™ Seal onto any crimp-top vial. It's that easy! The Poly Crimp™ Seal also forms an airtight seal by hand when the additional sealing of a crimping tool is not required. The Poly Crimp™ Seal offers the following benefits:

- Works on any crimp style vial.
- Crimps smoothly and easily—no more sharp edges from an aluminum crimp cap.
- Easily decaps using standard decapping tools.
- Eliminates inconsistent sealing that results from improper crimping.
- Unique design holds the septum firmly in place.
- Provides better seal than snap caps.*
- Meets CONEG standards. The Poly Crimp™ Seal is free of heavy metals and is easily disposed of or incinerated.



Description	100-pk.	400-pk.	1000-pk.
Clear Seal, PTFE/Red Rubber Liner	24433	—	24434
Clear Seal, PTFE/Silicone Liner	24443	—	24444
Blue Seal, PTFE/Silicone Liner	24445	—	24446
Clear Seal, PTFE/Silicone/PTFE Liner	24453	—	24454
Blue Seal, PTFE/Silicone/PTFE Liner	24455	—	24456
Mixed Color Seals, PTFE/Silicone/PTFE Liner	—	21728	—

*Based on continual use for 3–4 days.

tech tip

Compatibility concerns?

Refer to the Septum Selection Guide on [page 263](#).

2.0mL, 11mm GC Snap-Top Cap™ with Septum

- Clear polypropylene cap provides a tighter seal than standard Snap-Top Cap™.
- Good for pressurized injection systems.
- Save money! Use this snap cap when you don't need to crimp.



Description	100-pk.	1000-pk.
Polypropylene Cap, 10mil PTFE Liner	21729	21730
Polypropylene Cap, PTFE/Natural Rubber Septum	21731	21732
Polypropylene Cap, PTFE/Silicone Septum	21733	21734