

Convenience Kits: Vials, Caps, & Septa

Vials packaged in a clear-lid tray. Caps with septa packaged in a plastic bag.



Description	100-pk.	1000-pk.
2.0mL Clear Vial, deactivated , PTFE/Natural Rubber Silver Seal	24671	24672
2.0mL Amber Vial, deactivated , PTFE/Natural Rubber Silver Seal	24673	24674
2.0mL Clear Vial, untreated, PTFE/Natural Rubber Silver Seal	21196	21197
2.0mL Amber Vial, untreated, PTFE/Natural Rubber Silver Seal	21198	21199
2.0mL Clear Vial, untreated, PTFE/Silicone Silver Seal	24646	24647
2.0mL Amber Vial, untreated, PTFE/Silicone Silver Seal	24648	24649



USD List
24672, \$383.00
24674, 383
21197, \$255.00
21199, 255
24647, \$332.00
24649, 332

Limited Volume Inserts for 2.0mL Crimp-Top & Short-Cap, Screw-Thread Vials

Description	100-pk.	1000-pk.
50µL Glass, Big Mouth w/ Bottom Spring	24513	21782
250µL Glass, Big Mouth Insert w/ Bottom Spring	21776	21777
250µL Glass, Big Mouth Insert w/ Glass Flange (Step™ Design)*	24516	21779
350µL Glass, Flat Bottom Insert	21780	24517
350µL Glass, Flat Bottom Insert w/ ID Ring	24692	24693
250µL Polypropylene, Big Mouth, Bottom Spring	24518	—
250µL Polypropylene, Big Mouth, Top Flange	24519	—
250µL Polypropylene, Big Mouth, No Spring	24520	—

*Not to be used with 9mm screw-thread vials.



Limited Volume Vials, 11mm Crimp-Top

Description	100-pk.	1000-pk.
100µL Polypropylene	24651	24652
100µL Glass Insert/Clear Plastic Vial	24653	—



Versa Vial™ Vials, Glass, 2.0mL, 12 x 32mm

- Increased neck length for accommodating fluctuations in robotic arm settings.
- Eliminate cracked necks, breakage and dangerous sharp edges, using soft Versa Vial™ PTFE/silicone plugs.
- Largest target area available, ideal for sampling and pipet use.
- Vials available with or without graduated marking spots.
- Patent pending.

Description	100-pk.	1000-pk.
2.0mL Clear Glass Vial	22701	22702
2.0mL Clear Glass Vial w/Graduated Marking Spot	22703	22704
2.0mL Amber Glass Vial	22697	22698
2.0mL Amber Glass Vial w/Graduated Marking Spot	22699	22700



Versa Vial™ Vials, Polypropylene, 1.5mL, 12 x 32mm

Description	100-pk.	1000-pk.
1.5mL Polypropylene Vial	22705	22706



Versa Vial™ Inserts

Description	100-pk.	1000-pk.
250µL Glass Insert w/Glass Flange	22707	22708
250µL Polypropylene w/Flange	22709	22710
250µL Glass Flat Bottom Insert	22711	22712



Versa Vial™ Plugs

Description	100-pk.	1000-pk.
12mm, Green Polyethylene Plug	22713	22714
12mm, Grey Chlorobutyl/Siliconized Plug	22715	22716
12mm, White PTFE/Silicone Plug	22717	22718
12mm, White PTFE/Silicone Plug w/Slit*	22719	22720



VIALS & SYRINGES | VIALS