



Headspace Autosampler Vials

Description	100-pk.	1000-pk.
a) 6mL Clear Vial	21166	21167
b) 10mL Clear Vial, Flat Bottom	24683	24684
c) 10mL Clear Vial, Rounded Bottom	21164	21165
d) 20mL Clear Vial, Flat Bottom	24685	24686
e) 20mL Clear Vial, Rounded Bottom	21162	21163
f) 27mL Clear Vial	21160	21161

USD LIST	
21167	217
24684	266
21165	268
24686	280
21163	280
21161	408

a) - f) vial to instrument compatibility are designated in instrument reference chart below.

instrument reference

For Headspace Vials

The following instruments are compatible with the headspace vials listed above (as letter designated above). For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #	Manufacturer	Instrument/Model #
a) Carlo Erba	HS 500, HS 800	d) A.I.	Headspace
Agilent	Headspace	Carlo Erba	HS 500, HS 800
PerkinElmer	ISS 200	CTC	Headspace
		Cueni	Headspace
b) A.I.	Headspace	Dani	HSS 39.50, HSS 86.50, SPT 37.50
Carlo Erba	HS500, HS 800	Agilent	Headspace
Cueni	Headspace		
Dani	HSS 39.50, HSS 86.50	e) A.I.	Headspace
Agilent	Headspace	Cueni	Headspace
		Dani	HSS 39.50, HSS 86.50, SPT 37.50
c) A.I.	Headspace	L.E.A.P. Technologies	HS 500
Cueni	Headspace	PerkinElmer	Autosystem GC, F40, F45, HS 6, Integral 4000, HS 40
Dani	HSS 39.50, HSS 86.50	Tekmar	7000, 7000/7050
L.E.A.P. Technologies	HS 500	Varian	Genesis
PerkinElmer	Autosystem GC, F40, F45, HS 6, Integral 4000, HS 40		
Tekmar	7000, 7000/7050	f) Shimadzu	HSS-2B
Varian	Genesis		

tech tip

Compatibility concerns?
Refer to the Septum Selection Guide on page 263.

20mm Aluminum Seals w/Septa, Assembled



Description	100-pk.	1000-pk.
Silver Seal w/ PTFE/Gray Butyl Rubber	21761	21762
Silver Seal w/ PTFE/Silicone	21763	21764
Pressure Release Silver Seal w/ PTFE/Gray Butyl Rubber Septum	21765	21766
Pressure Release Silver Seal w/ PTFE/Silicone Septum	21767	21768
20mm Magnetic Crimp Cap w/Teflon®/Silicone Septum (tan/white) with 8mm Hole	22831	22832
20mm Magnetic Crimp Cap w/Teflon®/Silicone Septum (tan/white) with 5mm Hole	22833	22834
20mm Magnetic Crimp Cap w/PTFE/Butyl Septum with 8mm Hole	22835	22836
20mm Magnetic Crimp Cap w/PTFE/Butyl Septum with 5mm Hole	22837	22838

USD LIST	
21761	\$432.00
21763	\$369.00

5mm hole is compatible with the following systems: Carlo Erba HS500/HS800, CTC 500, Fisons HS500/HS800 and Varian/Chrompack 9020/25. 8mm hole is compatible with the Combi Pal.

also available

See page 274 for 20mm electronic crimper and decapper.

