

4.0mL Screw-Thread Step™ WISP™ 48 Style Vials, 13mm/425 Screw-Thread

4.0mL WISP™ 48 Screw-Thread Step™ Vials (15 x 45mm, 13/425 Screw-Thread)

Description	100-pk.	1000-pk.
4.0mL Clear Step™ Vial, White Grad. Marking Spot	24654	24655
4.0mL Amber Step™ Vial, White Grad. Marking Spot	24656	24657

Caps and Septa

Description	100-pk.	1000-pk.
Top Hat™, PTFE/Silicone Septa	21745	21746
Pre-Assembled, Red PTFE/Silicone, 0.065"	21743	21744
Caps, Black, Polypropylene, Open-Hole	21741	21742
Septa, Red PTFE/Silicone, 0.040"	21751	21752
Septa, Red PTFE/Silicone, 0.060"	21749	21750



4.0mL Clear Step™ Vial

Limited Volume Inserts

Description	100-pk.	1000-pk.
500µL, Glass, WISP™, Flat Bottom**	21787	21788

Polypropylene inserts available on request (1000-packs only).

**Also fit 4.0mL Snap Seal™ Vials (page 272).

1.0mL WISP™ 96 Style Shell Vials, 8 x 40mm

1.0mL WISP™ 96 Style, 8 x 40mm Shell Vials

Description	200-pk.	1000-pk.
1.0mL Clear Vial	24662	24663
1.0mL Amber Vial	24664	24665
1.0mL Polypropylene Vial	24678	24679
700µL Limited Volume Polypropylene Vial	24680	24681



1.0mL WISP™ 96 Style Vials

8mm Starburst Plugs

Description	200-pk.	1000-pk.
Snap Plug, Clear Polyethylene	21769	21770
Snap Plug, Conical, Clear Polyethylene	21771	21772
Snap Plug, Silicone-Lined, Clear Polyethylene	21773	21774

Limited Volume Inserts for Shell Vials

Description	100-pk.	1000-pk.
250µL Glass, Big Mouth Insert w/ Bottom Spring	21776	21777



Limited Volume Insert for Shell Vials

instrument reference

For 1.0mL WISP™ 96 Style Vials

The following instruments are compatible with 1.0mL WISP™ 96 style vials. For further information or questions, please contact your Restek sales representative or technical service.

Manufacturer	Instrument/Model #
Hitachi	L-7200
Waters	96-position M-700, 710, 717+

WHATMAN SYRINGE FILTERS

Whatman Syringe Filters Shipped Worldwide



WHATMAN GD/X syringe filters the best choice for the hard-to-filter samples. The WHATMAN syringe filters are made with a pigment-free polypropylene housing to eliminate sample contamination. WHATMAN syringe filters contain a unique prefiltration stack of Whatman GMF and GF/F glass microfibre media allowing filtration of difficult samples with less hand pressure.

Whatman Syringe Filters can be used in conjunction with the Cornwall® Syringe fluid dispensing system below.



- * Sterile- all materials used in constructions pass USP Class VI
- * 25mm filters below are for samples equal or greater than 10mL
- * Individually blister packed
- * Sold as a pack of 10 or 50

[CLICK HERE TO ORDER](#)

Item No.	Description	Sale Unit	Price USD
6900-2504-50	Whatman 25 mm GD/X Syringe Filters Contain a unique prefiltration stack of Whatman multigrade GMF 150 (10:1 um) and Grade GF/F (0.7 um) prefilters. Final filter has a 0.45 um filter pore and PVDF (PolyVinylidene Fluoride) membrane	Pack of 50	\$142.00
		Pack of 10	\$40.00
6902-2504-50	Similar to the above yet without the GF/F prefilter. Final filtering membrane is glass microfiber	Pack of 50	\$142.00
		Pack of 10	\$40.00

CORNWALL® FLUID DISPENSING SYRINGE



The Cornwall® Syringe fluid dispensing system is effectively used when repetitious filling of a syringe to the same volume is needed. It is a manual, sterile device that has a two-way valve that attaches to both the syringe and the solution to be drawn into the syringe. Each time the grip is released, the syringe fills with solution. Each time the grip is compressed, the syringe volume is expelled. The Cornwall Syringe System is a sterile, disposable system that includes spring-loaded syringe, valve system, spike and tubing set for repetitive fluid dispensing. 10cc capacity.

CAN BE USED IN CONJUNCTION WITH WHATMAN SYRINGE FILTERS ABOVE.

Cornwall Syringe Fluid Dispensing System

10mL, Sterile, Disposable

Item No. 305224

Price Each: \$20.00 USD



VIAL SEAL
Center Tear Out
15 colors & stripes
13mm, 20mm & 30mm



VIAL SEAL
Aluminum Tear Off
9 colors & stripes
13mm, 20mm & 30mm



VIAL SEAL
Aluminum with Flip Cap
17 colors
13mm & 20mm



VIAL SEAL
Tear Off Flip Cap
4 colors
13mm & 20mm



SYRINGE FILTERS



BOTTLE FILTERS



VIAL STOPPERS
13mm and 20mm
Round Bottom & 2-Leg



VIAL CRIMPERS
For all 13mm and 20mm
Vial Seals



SERUM VIALS
Clear and Amber
2mL thru 200mL

VOIGT GLOBAL DISTRIBUTION INC

PO Box 1130
Lawrence, Kansas 66044 USA
Toll Free in the USA: 877-484-3552 or 1.785.749.6777
info@vial-seals.com