

Electronic Crimpers and Decappers

Hand-held Rechargeable Crimpers and Decappers

- Easy to use; comfortable grip.
- Hundreds of operations from one charge.
- Adjustable crimping force.
- New battery allows for faster charge cycle.
- All kits include a universal plug set for the battery charger.



new and improved!

All electronic crimpers & decappers now come with a universal plug set, which includes U.S., European, UK, and Australian plug clips.

The newly re-designed electronic crimpers and decappers have a 4.8 volt battery already installed. The new battery allows for quick, powerful operation and a faster charge cycle, minimizing downtime. Charging the battery now only takes thirty minutes to two hours. Also, the power supply comes with a plug set for operation in most countries. The universal plug set includes U.S., European, UK and Australian plug clips.

These hand-held electronic tools for the 11mm and 20mm caps will cycle hundreds of times on a single charge. The cycle is controlled precisely, through an internal counter, and is adjusted with two buttons on the side of the case. The crimpers and decappers easily fit around vials in a standard autosampler tray. Each kit includes a tool, a 4.8v battery, which allows for a faster charge cycle, and a universal plug set for the battery charger.



Capping
Easily crimp caps on standard crimp vials.

Decapping has never been easier!

one...



two...



three!



Universal plug set for battery charger included in all kits.

Description	qty.	cat.#
11mm Electronic Crimper	kit	22338
11mm Electronic Decapper	kit	22340
20mm Electronic Crimper	kit	22339
20mm Electronic Decapper	kit	22341
Replacement Battery for Electronic Crimpers and Decappers	ea.	22342