

Aluminum Seal Crimpers and Decappers

The crimper is adjustable for optimized sealing performance and is comfortable enough for prolonged use. If you need to save, transfer, or dispose of your samples, our decapper allows you to remove a crimp-top cap safely and easily. Far safer than removing caps with pliers, knife or screwdriver blades, etc. If you aren't currently using a seal decapper, we recommend you order one today!

Description	qty.	cat.#
8mm Crimper	ea.	21735
8mm Decapper	ea.	21736
11mm Crimper	ea.	21170
11mm Decapper	ea.	21171
13mm Crimper	ea.	21739
13mm Decapper	ea.	21740
20mm Crimper	ea.	21737
20mm Decapper	ea.	21738



Safer than removing caps with pliers, knives or screwdriver blades!



Decapper



Crimper

Snap Rack™ Storage Unit

The Snap Rack™ sample vial rack is an innovation that allows you to keep your vials and caps together. Constructed of tough polypropylene, the Snap Rack™ holds 50 12 x 32mm size vials. The Snap Rack™ is stackable too! Simply “snap” one rack on top of another. Double rack capacity by snapping a rack across another rack. With a standard vial rack, once the post breaks, a replacement rack is needed. The Snap Rack™ eliminates this problem!



Snap Rack™ vial racks are stackable for maximum storage of vials and caps!

Description	Capacity	qty.	cat.#
50 Position Snap Rack™ for 12 x 32mm Vials	50	ea.	24521

Mininert® Precision Sampling Valves for Vials

Mininert® valves are on/off valves that are very convenient for repetitive sampling, while limiting sample exposure to the silicon septum. Models are available for screw-cap and crimp-top vials. The screw-cap Mininert® valve is available in 13, 15, 18, 20, and 24mm cap sizes. The crimp-top valve for 20mm ID glassware slides into the neck of the vial. Turn the threaded flange to secure a tight fit.

Description	Thread Size	qty.	cat.#
13mm Mininert® Screw Cap	13mm-425	12-pk.	24900
15mm Mininert® Screw Cap	15mm-425	12-pk.	24901
18mm Mininert® Screw Cap	18mm-400	12-pk.	24902
20mm Mininert® Screw Cap	20mm-400	12-pk.	24903
24mm Mininert® Screw Cap	24mm-400	12-pk.	24904
20mm Mininert® Crimp Top	—	12-pk.	24905
Replacement Septa	—	50-pk.	24906
Septum Insertion Tool	—	ea.	24907



Mininert® Screw Cap



Mininert® Crimp Top

also available

See page 269 for Accuform Micro-Vials to use with Mininert® caps & valves.

