

Autosampler Syringes for Varian GCs

Autosampler Syringes for Varian 8000 Series GCs



Steve Kozel

U.S. Sales & Distribution
Manager
1+ year of service!



Hamilton 10µL Autosampler Syringe

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton Model	Hamilton cat.#	qty.	Restek cat.#
Varian 8100/8200 Syringes								
10µL	RN	26s	1.99"	5	701 Varian	202880	ea.	24956



SGE 10µL Autosampler Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE Model	SGE cat.#	qty.	Restek cat.#
Varian 8000 Syringes								
10µL	R	25	50mm	Dome	10R-VA80-II	002926	ea.	24850
Varian 8035, 8100, & 8200 Syringes								
10µL*	F	25	53mm	S/Hole	10F-VA8X-11	002923	ea.	24851
10µL*	R	25	53mm	S/Hole	10R-VA8X-II	002924	ea.	24852
100µL	R	25	53mm	S/Hole	100R-VA8X	005921	ea.	24853
Varian CP 8400, CP8410 Syringes								
10µL	F	26	50mm	Cone	10F-C/F-5/0.47C	002980	ea.	24922
10µL	F	26	50mm	Bevel	10F-VA8400-5/0.47	002950	ea.	21202

*Gas-tight syringe.

Replacement Needles for Varian 8000 Series Autosampler Syringes

Needles for SGE Syringes

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
10µL	23	53mm	S/Hole	037779	2-pk.	24857
10µL*	25	53mm	S/Hole	037780	ea.	24858
10µL**	25	53mm	S/Hole	037777	ea.	24854
10µL	25	53mm	2	037776	2-pk.	24855

*straight needle shaft, 0.20mm ID needle, for use with viscous samples

**tapered needle shaft, 0.12mm ID needle

Searching for a product?

www.restek.com