

Standard Microliter Volume Syringes

Standard Microliter Volume Syringes

- Economical.
- Available with cemented needles (N/F) or removable needles (RN/R).
- Each syringe plunger and barrel assembly is manufactured as one working unit; components are not interchangeable or individually replaceable.



Hamilton 700 Series Microliter Syringe

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		qty.	Restek cat.#
					Model	cat.#		
5 μ L	N	26s	2"/51mm	2	75N	87900	ea.	24938
5 μ L	RN	26s	2"/51mm	2	75RN	87930	ea.	24617
10 μ L	N	26s	2"/51mm	2	701N	80300	ea.	20174
10 μ L	N	26s	2"/51mm	2	701N	80366	6-pk.	20175
10 μ L	RN	26s	2"/51mm	2	701RN	80330	ea.	24530
25 μ L	N	22s	2"/51mm	2	702N	80400	ea.	24531
10 μ L	N	26s	2"/51mm	5	701NPT5	80339	ea.	24967
25 μ L*	RN	22s	2"/51mm	2	702RN	80430	ea.	24532
50 μ L	N	22s	2"/51mm	2	705N	80500	ea.	24533
50 μ L*	RN	22s	2"/51mm	2	705RN	80530	ea.	24534
100 μ L	N	22s	2"/51mm	2	710N	80600	ea.	24535
100 μ L*	RN	22s	2"/51mm	2	710RN	80630	ea.	24536
250 μ L	N	22s	2"/51mm	2	725N	80700	ea.	24537
250 μ L*	RN	22s	2"/51mm	2	725RN	80730	ea.	24538
500 μ L	N	22	2"/51mm	2	750N	80800	ea.	24539
500 μ L*	RN	22	2"/51mm	2	750RN	80830	ea.	24540



SGE Standard Microliter Volume Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		qty.	Restek cat.#
					Model	cat.#		
5 μ L	F	26	50mm	2	5F**	001000	ea.	24700
5 μ L*	R	26	50mm	2	5R**	001050	ea.	24701
10 μ L	F	26	50mm	2	10F**	002000	ea.	24702
10 μ L	F	26	50mm	2	SK-10F**	002030	6-pk.	24715
10 μ L*	R	26	50mm	2	10R**	002050	ea.	24703
10 μ L*	R	26	50mm	2	SK-10R**	002080	6-pk.	24716
25 μ L	F	25	50mm	2	25F	003000	ea.	24704
25 μ L*	R	25	50mm	2	25R	003050	ea.	24705
50 μ L	F	25	50mm	2	50F	004000	ea.	24706
50 μ L*	R	25	50mm	2	50R	004050	ea.	24707
100 μ L	F	25	50mm	2	100F	005000	ea.	24708
100 μ L*	R	25	50mm	2	100R	005050	ea.	24709
250 μ L	F	25	50mm	2	250F	006000	ea.	24710
250 μ L*	R	25	50mm	2	250R	006050	ea.	24711
500 μ L	F	25	50mm	2	500F	007000	ea.	24712
500 μ L*	R	25	50mm	2	500R	007050	ea.	24713

*For replacement RN/R needles, see page 290.

**With plunger protection.

