

SuperfleX Flexible Plunger Syringes

- The SuperfleX syringe features a super-elastic plunger that returns to its original shape if twisted or bent.



SGE SuperfleX Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE			Restek cat.#
					Model	cat.#	qty.	
5 μ L	F	26	50mm	2	5FX	001100	ea.	24720
5 μ L	R	26	50mm	2	5RX	001150	ea.	24721
10 μ L	F	26	50mm	2	10FX	002100	ea.	24722
10 μ L	F	26	50mm	2	SK-10FX	002130	6-pk.	24718
10 μ L	R	26	50mm	2	10RX	002150	ea.	24723
10 μ L*	R	26	50mm	2	SK-10RX	002180	6-pk.	24719

Reinforced Plunger Microliter Volume Syringes

- Reinforcing extends the life of fragile, small diameter plungers.
- Economical because of longer life—even in rugged applications.
- Hand-lapped plunger ensures a perfect seal with the barrel.



Hamilton 900 Series Microliter Syringe

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton			Restek cat.#
					Model	cat.#	qty.	
5 μ L	N	26s	2"/51mm	2	95	87920	ea.	24541
5 μ L*	RN	26s	2"/51mm	2	95	87925	ea.	24542
10 μ L	N	26s	2"/51mm	2	901	80360	ea.	24543
10 μ L*	RN	26s	2"/51mm	2	901	80370	ea.	24544

Guide to Needle Termination Codes

Hamilton:

(N) Cemented Needle
 (RN) Removable Needle
 (ASN) Autosampler Cemented Needle
 (ASRN) Autosampler Removable Needle
 (TLL) Teflon® Luer Lock
 (KH) Knurled Hub
 (LT) Luer Tip
 (LTN) Luer Tip Cemented Needle

SGE:

(F) Fixed Needle
 (R) Removable Needle
 (LL) Luer Lock
 (LT) Luer Tip



SGE Guided Plunger Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE			Restek cat.#
					Model	cat.#	qty.	
5 μ L	F	26	50mm	2	5F-GP	001400	ea.	24769
5 μ L*	R	26	50mm	2	5R-GP	001450	ea.	24766
10 μ L	F	26	50mm	2	10F-GP	002400	ea.	24767
10 μ L*	R	26	50mm	2	10R-GP	002450	ea.	24768

*For replacement RN/R needles, see page 290.