

Teflon® Tip, Gas-Tight Syringes

- Suitable for gases or liquids for maximum inertness.
- High accuracy of dispensed volumes.
- Interchangeable barrels, plungers, and tips extend performance and increase cost-effectiveness.



Hamilton 1700/1000 Series Gas-Tight Syringe

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek
					Model	cat.#	
10µL	N	26s	2"/51mm	2	1701	80000	ea. 24557
10µL	RN	26s	2"/51mm	2	1701	80030	ea. 24558
25µL	N	22s	2"/51mm	2	1702	80200	ea. 24559
25µL	RN	22s	2"/51mm	2	1702	80230	ea. 24560
50µL	N	22s	2"/51mm	2	1705	80900	ea. 24561
50µL	RN	22s	2"/51mm	2	1705	80930	ea. 24562
100µL	N	22s	2"/51mm	2	1710	81000	ea. 24563
100µL	RN	22s	2"/51mm	2	1710	81030	ea. 24564
250µL	N	22s	2"/51mm	2	1725	81100	ea. 24567
250µL	RN	22s	2"/51mm	2	1725	81130	ea. 24568
500µL	LTN	22	2"/51mm	2	1750	81217	ea. 24571
500µL	RN	22	2"/51mm	2	1750	81230	ea. 24572
1mL	LTN	22	2"/51mm	2	1001	81317	ea. 24575
1mL	RN	22	2"/51mm	2	1001	81330	ea. 24576
1mL	TLL	*	—	w/out slots	1001	81320	ea. 24578
2.5mL	TLL	*	—	w/out slots	1002	81420	ea. 24584
2.5mL	RN	22	2"/51mm	2	1002	81430	ea. 24582
2.5mL	N	22	2"/51mm	2	1002	81417	ea. 24581
5mL	TLL	*	—	w/out slots	1005	81520	ea. 20178
10mL	TLL	*	—	w/out slots	1010	81620	ea. 20179
25mL	TLL	*	—	w/out slots	1025	82520	ea. 24591
50mL	TLL	*	—	w/out slots	1050	85020	ea. 24717



SGE Gas-Tight Syringe

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek
					Model	cat.#	
10µL	F	26	50mm	2	10F-GT	002200	ea. 24725
10µL	R	26	50mm	2	10R-GT	002250	ea. 24726
25µL	F	25	50mm	2	25F-GT	003200	ea. 24727
25µL	R	25	50mm	2	25R-GT	003250	ea. 24728
50µL	F	25	50mm	2	50F-GT	004200	ea. 24729
50µL	R	25	50mm	2	50R-GT	004250	ea. 24730
100µL	F	25	50mm	2	100F-GT	005200	ea. 24734
100µL	R	25	50mm	2	100R-GT	005250	ea. 24735
100µL	LL	*	—	—	100F-LL-GT	005230	ea. 24737
250µL	F	25	50mm	2	250F-GT	006200	ea. 24738
250µL	R	25	50mm	2	250R-GT	006250	ea. 24739
250µL	LL	*	—	—	250F-LL-GT	006230	ea. 24741
500µL	F	25	50mm	2	500F-GT	007200	ea. 24742
500µL	R	25	50mm	2	500R-GT	007250	ea. 24743
1mL	R	23	50mm	2	1MR-GT	008100	ea. 24750
1mL	LL	*	—	—	1MDF-LL-GT	008025	ea. 24752
2.5mL	LL	*	—	—	2.5MDF-LL-GT	008425	ea. 24755
5mL	LL	*	—	—	5MDR-LL-GT	008760	ea. 24757
10mL	LL	*	—	—	10MDR-LL-GT	008960	ea. 24759
25mL	LL	*	—	—	25MDR-LL-GT	009462	ea. 24760
50mL	LL	*	—	—	50MR-LL-GT	009660	ea. 24761
100mL	LL	*	—	—	100MR-LL-GT	009760	ea. 24762

*Needles sold separately. See page 290.

Guide to Needle Termination Codes

Hamilton:
 (N) Cemented Needle
 (RN) Removable Needle
 (ASN) Autosampler
 Cemented Needle
 (ASRN) Autosampler
 Removable Needle
 (TLL) Teflon® Luer Lock
 (KH) Knurled Hub
 (LT) Luer Tip
 (LTN) Luer Tip Cemented
 Needle

SGE:
 (F) Fixed Needle
 (R) Removable Needle
 (LL) Luer Lock
 (LT) Luer Tip

