



Alan Sensus
Technical Service
Specialist
6+ years of service!

A-2 Luer Gas-Tight Syringes*

- For use with VOA-CYL 4 standards.
- Push-style on-off button valve.

Use to withdraw a gas sample from a high-pressure cylinder through a high-purity VOC single-stage regulator after attaching the appropriate adaptor to the regulator. Use with luer syringe needles listed.



100µL Luer Syringe

Description	qty.	cat.#
100µL Luer Syringe	ea.	20162
500µL Luer Syringe	ea.	20163
Luer Syringe Needles (26s gauge)	6-pk.	20133

*Restek's A-2 syringes have a luer tip.

Heavy Duty Purge & Trap Syringe (Dynatech Precision Sampling)

- Heavy-duty glass barrel with metal front and rear flanges.
- Teflon® luer-lock tip.
- Can fill and empty sparge tubes.
- Sample-Lok valve accepts standard luer-lock needles.
- Suitable for purge & trap applications.



Precision Sampling Syringe
with Sample-Lok Valve

Syringe	cat.#	2mL	5mL	10mL
without Sample-Lok	21205		21206	21209
with Sample-Lok	—		21208	21207

Luer-Lock Valve Adaptors



Description	SGE cat. #	qty.	Restek cat.#
Luer-Lock Valve Adaptors	200010	2-pk.	21283

Jumbo Syringe

Clear acrylic syringes, ideal for holding and dispensing large volumes of gas. An adjustable plunger on the O-ring ensures that the syringe is gas-tight over a long period of time. The central port is supplied with a luer-lock fitting; the secondary port is supplied with a septum nut. This enables access to the gas sample for adding standards or removing a sub-sample. The plunger stem is detachable, making sample storage easy.



Jumbo Syringe

1000mL Jumbo Syringe



O-Rings for 1000mL & 500mL Syringes

Volume	Model	SGE cat. #	qty.	Restek cat.#
500mL	500MAR-LL-GT	009910	ea.	21275
1000mL	1000MAR-LL-GT	009920	ea.	21276
2000mL	2000MAR-LL-GT	009930	ea.	21277

Syringe O-Rings

Syringe Volume	SGE cat.#	qty.	Restek cat.#
500mL	032527	ea.	21278
1000mL	032532	ea.	21279