

**Matt Matis**

Northeast States Sales
Representative
4+ years of service!

Gas-Tight Syringes for Rheodyne® & Valco® Valves

- Teflon®-tipped plungers.
- Removable needles.
- Replaceable syringe barrels, plungers, and plunger tips.



Hamilton 10µL Gas-Tight Syringe for
Rheodyne® and Valco® Valves

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		Restek
					Model	cat.#	cat.#
10µL	RN	26s	2"/51mm	3	1701	80065	ea. 21260
25µL	RN	22s	2"/51mm	3	1702	80265	ea. 21261
50µL	RN	22s	2"/51mm	3	1705	80965	ea. 21262
100µL	RN	22s	2"/51mm	3	1710	81065	ea. 21263
250µL	RN	22s	2"/51mm	3	1725	81165	ea. 21264



SGE 50µL Gas-Tight Syringe for Rheodyne® and Valco® Valves

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		Restek
					Model	cat.#	cat.#
10µL	R	22	2"	3	10R-GT-LC	002313	ea. 24866
25µL	R	22	2"	3	25R-GT-LC	003312	ea. 24867
50µL	R	22	2"	3	50R-GT-LC	004312	ea. 24868
100µL	R	22	2"	3	100R-GT-LC	005312	ea. 24869
250µL	R	22	2"	3	250R-GT-LC	006312	ea. 24870
500µL	R	22	2"	3	500R-GT-LC	007312	ea. 24871

Replacement Needles for Gas-Tight Syringes for Rheodyne® & Valco® Valves



Needles for Hamilton Syringes

Syringe Volume	Needle Gauge	Needle Length	Point Style	Hamilton cat.#	qty.	Restek cat.#
5–100µL	22s	2"	3	7770-01	6-pk.	24941
250µL–10mL	22	2"	3	7780-04	6-pk.	24945

Needles for SGE Syringes

Syringe Volume	Needle Gauge	Needle Length	Point Style	SGE cat.#	qty.	Restek cat.#
10µL	22	2"	3	037250	5-pk.	24808
25–500µL	22	2"	3	038250	5-pk.	24809