

Economical Microliter Liquid-Type Syringes for Rheodyne® & Valco® Valves

- Cemented/fixed needles.
- Plungers and barrels are not interchangeable or replaceable.



Hamilton 50µL Microliter Syringe for Rheodyne® and Valco® Valves

Hamilton Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	Hamilton		qty.	Restek cat.#
					Model	cat.#		
10µL	N	22s	2"	3	701	80365	ea.	21250
25µL	N	22s	2"	3	702	80465	ea.	21251
50µL	N	22s	2"	3	705	80565	ea.	21252
100µL	N	22s	2"	3	710	80665	ea.	21253
250µL	N	22s	2"	3	725	80765	ea.	21254



SGE 50µL Microliter Syringe for Rheodyne® and Valco® Valves

SGE Syringes

Volume	Needle Term.	Needle Gauge	Needle Length	Point Style	SGE		qty.	Restek cat.#
					Model	cat.#		
10µL	F	22	2"	3	10F-LC	002301	ea.	24860
25µL	F	22	2"	3	25F-LC	003300	ea.	24861
50µL	F	22	2"	3	50F-LC	004300	ea.	24862
100µL	F	22	2"	3	100F-LC	005300	ea.	24863
250µL	F	22	2"	3	250F-LC	006300	ea.	24864
500µL	F	22	2"	3	500F-LC	007300	ea.	24865

Guide to Needle Termination Codes

Hamilton:

- (N) Cemented Needle
- (RN) Removable Needle
- (ASN) Autosampler Cemented Needle
- (ASRN) Autosampler Removable Needle
- (TLL) Teflon® Luer Lock
- (KH) Knurled Hub
- (LT) Luer Tip
- (LTN) Luer Tip Cemented Needle

SGE:

- (F) Fixed Needle
- (R) Removable Needle
- (LL) Luer Lock
- (LT) Luer Tip

Priming Syringes for Waters 6000, 6000A, 501, 510, 590, 610 and 610E HPLC Pumps

- Designed for maximum safety with metal flange and luer-lock (LL) hub.
- Teflon®-tipped plunger.



Hamilton Syringe

Volume	Needle Term.	Hamilton Model	Hamilton cat.#	Pressure Tested to	qty.	Restek cat.#
10mL	TLL	1010W	81610	700 psig	ea.	21265

SGE Syringe

Volume	Needle Term.	SGE Model	SGE cat.#	Pressure Tested to	qty.	Restek cat.#
10mL	LL	10MDR-LL-GT	008960	100 psig	ea.	24759