



### Bulk Packing Materials

#### Use our bulk packing materials to pack your own columns!

Restek is among the small group of column manufacturers capable of producing their own high purity silica. We offer our Pinnacle™ II and Pinnacle™ DB silica and bonded phases in bulk. Pinnacle™ II is an excellent replacement for Hypersil® silica. Pinnacle™ DB is a highly base-deactivated material for analyses of basic compounds and an excellent alternative to Hypersil® BDS silica. Bonded phases from our Ultra, Allure®, and Viva lines also are available in bulk.

Our extensive QC program ensures the high quality and reproducibility of these silicas. Each lot of material is tested for mean particle size and distribution, pore diameter, surface area, and total metals analysis. You can be confident that you are getting consistent, high-quality product.

Use these materials for easy scale up to preparative scale chromatography, or for packing your own columns.

### Bulk Packing Materials

Description	min. qty.	cat.#
<b>3µm Pinnacle™ II Bulk Packing Materials</b>		
Pinnacle™ II C8 Bulk Packing	5g	92133
Pinnacle™ II C18 Bulk Packing	5g	92143
Pinnacle™ II Cyano Bulk Packing	5g	92163
Pinnacle™ II Phenyl Bulk Packing	5g	92153
<b>5µm Pinnacle™ II Bulk Packing Materials</b>		
Pinnacle™ II Amino Bulk Packing	5g	92175
Pinnacle™ II C8 Bulk Packing	5g	92135
Pinnacle™ II C18 Bulk Packing	5g	92145
Pinnacle™ II Cyano Bulk Packing	5g	92165
Pinnacle™ II Phenyl Bulk Packing	5g	92155
<b>5µm Pinnacle™ DB Bulk Packing Materials</b>		
Pinnacle™ DB C18 Bulk Packing	5g	94145
Pinnacle™ DB C8 Bulk Packing	5g	94135
Pinnacle™ DB Cyano Bulk Packing	5g	94165
Description	qty.	cat.#
<b>5µm Ultra Bulk Packing Materials</b>		
Ultra C1 Bulk Packing	10g/btl.	91015
Ultra C4 Bulk Packing	10g/btl.	91025
Ultra C8 Bulk Packing	10g/btl.	91035
Ultra C18 Bulk Packing	10g/btl.	91745
Ultra Amino Bulk Packing	10g/btl.	91075
Ultra Cyano Bulk Packing	10g/btl.	91065
Ultra Phenyl Bulk Packing	10g/btl.	91055

### Silica Bulk Materials

Description	min. qty.	cat.#
<b>Silica Bulk Materials</b>		
Pinnacle™ II Silica Bulk Packing, 3µm	5g	92103
Pinnacle™ II Silica Bulk Packing, 5µm	5g	92105
Pinnacle™ DB Silica Bulk Packing, 5µm	5g	94105
Description	qty.	cat.#
Ultra Silica Bulk Packing, 5µm	10g/btl.	91005

### Kromasil® Bulk Packings

- High-purity packing materials in 10 and 16µm.
- All Kromasil® phases available, please inquire.

Description	min. qty.	cat.#
Kromasil® 100Å Silica, 10µm	200g	92000
Kromasil® 100Å Silica, 16µm	200g	92001
Kromasil® 100Å C8, 10µm	200g	92030
Kromasil® 100Å C8, 16µm	200g	92031
Kromasil® 100Å C18, 10µm	200g	92040
Kromasil® 100Å C18, 16µm	200g	92041
Kromasil® 100Å Chiral DMB, 10µm	200g	92080
Kromasil® 100Å Chiral DMB, 16µm	200g	92081
Kromasil® 100Å Chiral TBB, 10µm	200g	91990
Kromasil® 100Å Chiral TBB, 16µm	200g	91991

