

## Restek Replacement Parts for Agilent HPLC Systems



Piston Seals



Outlet Ball Valve



Needle Seat Assembly



Lamp, VWD G1314A

new!

new!

Description	Model #	Similar to Agilent part #	qty.	cat.#
Preventive Maintenance Kit (Includes: rotor seal, needle seat, needle assembly, seat cap)	1050	01078-68721	kit	25259
Autosampler Preventive Maintenance Kit (Includes: rotor seal, needle assembly, needle seat)	1100	G1313-68709	kit	25271
Pump Maintenance Kit (Includes: PTFE frit, outlet cap, active inlet cartridge, gold disk seal, 2 piston seals, glass solvent filter)	1050, 1100	G1311-68710	kit	25270
Outlet Ball Valve, Binary Pump	1100, 1200	G1312-60012	ea.	25267
Outlet Ball Valve	1050, 1100, 1200	G1311-60012	ea.	25276
Sieves for Outlet Valve	1050, 1100, 1200	5063-6505	10-pk.	25266
Check Valve Cartridge Assembly	1090	79835-67101	ea.	25344
Piston Seals, Teflon® w/Graphite	1050, 1100, 1200	5063-6589	2-pk.	22482
Piston Seals, Teflon® w/Graphite	1050, 1100, 1200	5063-6589	10-pk.	22483
Piston Seals (Black)	1090	5062-2494	4-pk.	25347
Seal Wash Kit, Binary Pump (4 seals, 4 gaskets)	1100	—	kit	25268
Seal Wash Kit (2 seals, 2 gaskets)	1050, 1100	—	kit	25269
Wash Seal	1050, 1100, 1200	0905-1175	ea.	25277
Sapphire Piston	1050, 1100, 1200	5063-6586	ea.	25273
Sapphire Piston	1090	6980-0672	ea.	25345
Needle Seat	1050	79846-67101	ea.	25258
Needle Seat	1090	79846-67101	ea.	25348
Needle Seat Assembly	1100	G1313-87101	ea.	25265
Needle Assembly	1100	G1313-87201	ea.	25278
Rotor Seal (not for use with 7125 injection valve)	1050	0101-0626	ea.	25272
Rotor Seal	1100	0100-1853	ea.	25275
Rotor Seal (Rheodyne®-style)	1090	0101-0623	ea.	25349
Frits, PTFE	1050, 1100, 1200	01018-22707	5-pk.	25466
Seal, Gold Disk (outlet)	1050, 1100, 1200	5001-3707	ea.	25467
Outlet Cap	1050, 1100	5062-2485	4-pk.	25139
Outlet Cap & Gold Seal Assembly	1050, 1100	—	2-pk.	25140
Connecting Tube	1050, 1100	G1311-67304	ea.	25058
Bottle Head Assembly	—	G1311-60003	ea.	25059
Detector Lamp, 1090 DA, 1050 VW/DA/MWD	1090, 1050	79883-60002	ea.	25260
Lamp, DAD G1315A, G1365A	1100, 1200	2140-0590	ea.	25261
Lamp, VWD G1314A	1100, 1200	G1314-60100	ea.	25262
8453 Deuterium Lamp	—	2140-0605	ea.	25263
G1321 Fluorescence Detector Flash Lamp	—	2140-0600	ea.	25264
Lamp, DAD Long Life Deuterium (2000 hours)	1100	5181-1530	ea.	25399
Helium Tubing, 5m, PTFE, 0.057" ID	1050, 1100	5062-2461	ea.	24990
Injector Lubricant, PTFE, 10g	1050, 1100	79841-65501	ea.	26520

## Capillary Stainless Steel Tubing Assemblies for Agilent HPLC Systems



Capillary SS Tubing with fittings



Mixing Capillary Assembly

Description	Model #	Similar to Agilent part #	qty.	cat.#
Capillary SS Tubing, 130mm x 0.17mm ID, with fittings	1090	01090-87305	ea.	26525
Capillary SS Tubing, 800mm x 0.17mm ID, with fittings	1050	01078-87305	ea.	26526
Capillary SS Tubing, 180mm x 0.17mm ID, with fittings	1100	G1313-87305	ea.	26527
Capillary SS Tubing, 700mm x 0.25mm ID, with fittings	1050	01018-67305	ea.	26528
Capillary SS Tubing, 700mm x 0.25mm ID, with fittings	1050	01078-87306	ea.	26529
Seat Capillary, SS Tubing, 0.17mm ID	1050	01078-87303	ea.	26530
Capillary SS Tubing, 105mm x 0.17mm ID	—	5021-1816	ea.	26531
Mixing Capillary Assembly	1100	G1312-67302	ea.	26532
Capillary SS Tubing, Valve to Metering Head	1100	G1313-87301	ea.	26533
Capillary SS Tubing, 150mm x 0.17mm ID	—	5021-1817	ea.	26534
Capillary SS Tubing, 280mm x 0.17mm ID	—	5021-1818	ea.	26535
Capillary SS Tubing, 400mm x 0.17mm ID	—	5021-1819	ea.	26536
1/16" SS Fitting, Front and Back Ferrules	—	5062-2418	10-pk.	26537