

**Opti-Cap™ Bottle Top**

The most economical way to helium-sparge and deliver HPLC mobile phases. The Opti-Cap™ top fits all standard GL45 bottles and has two 1/8-inch holes and one 1/16-inch hole for tubing.

Description	qty.	cat.#
Opti-Cap™ (Cap and PEEK™ Plug)	ea.	25300
Opti-Cap™ Kit (Opti-Cap™, 3 meters of tubing, sparging filters)	kit	25301
Opti-Cap™ Kit with 1L Bottle	kit	25302
Opti-Cap™ Kit with 2L Bottle	kit	25303
Related items and replacement parts	qty.	cat.#
SS Mobile Phase Mobile Phase Sparge Filter: 2µm	ea.	25311
Mobile Phase Inlet Filter: 10µm	ea.	25312
Teflon® Tubing, 1/8" OD x 0.094" ID x 3m (2.4mm ID)	3m	25307
Teflon® Tubing, 1/8" OD x 0.063" ID x 3m (1.6mm ID)	3m	25306
PEEK™ Plug, 1/4"-28 threads	3-pk.	25319
1L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25304
2L Graduated Safety-Coated Bottle – GL-45 threads	ea.	25305



Opti-Cap™ Kit with bottle



**Eco-Cap Bottle Top, PTFE**

Fits all standard GL-45 bottles; has two 1/8-inch holes and one 1/16-inch hole for tubing.

Description	qty.	cat.#
Eco-Cap Bottle Top, includes 1 Male Luer Plug	ea.	25395

also **available**

See **page 352** for Hub-Cap bottle tops, adapters, and filter kits.



**Technical Service Specialists**

Al Carusone, Tom Bloom, Jonathan Keim, Tim Herring, Alan Sensue, Terry Reid