



PEEK™ Survival Kit



**Terry Reid**  
Technical Service  
Specialist  
9+ years of service!

### PEEK™ Survival Kit for HPLC

The PEEK™ Survival Kit is an invaluable parts kit that contains tubing, fittings, and tools essential for setting up and maintaining your HPLC system: PEEK™ tubing, connectors, and elbows, Teflon® tubing, a tubing cutter and extra blades, a ValvTool wrench, open-end wrenches, and more.

#### For start-up and maintenance in all HPLC systems.

Description	qty.	cat.#
PEEK™ Survival Kit for HPLC	kit	25322

#### Contents:

Description	qty.	cat.#
PEEK™ Column Connector (beige, round body)	10-pk.	25015
PEEK™ Tubing, 1/16" OD x 0.005" ID x 3m (red stripe)	ea.	25065
PEEK™ Tubing, 1/16" OD x 0.007" ID x 3m (yellow stripe)	ea.	25066
PEEK™ Tubing, 1/16" OD x 0.010" ID x 3m (blue stripe)	ea.	25067
PEEK™ Tubing Elbow 90°	5-pk.	25308
PEEK™ Tubing Elbow 180°	5-pk.	25309
Teflon® Tubing, 1/8" OD x 0.063" ID x 3m (1.6mm ID)	3m	25306
Teflon® Tubing, 1/8" OD x 0.094" ID x 3m (2.4mm ID)	3m	25307
Tubing Clip	5-pk.	25310
ValvTool Wrench	ea.	25321
Open-End Wrenches (1/4" x 5/16")	2-pk.*	20110
Clean-Cut™ Tubing Cutter	ea.	25069
Replacement Blade for Clean-Cut™ Cutter	ea.	25070
PEEK™ Union Connector 1/16"	2-pk.	25323
Mobile Phase Sparge Filter: 2µm	ea.	25311
Inlet Filter: 10µm	ea.	25312

\*Kit contains 1 wrench, replacement (cat.# 20110) is a 2-pk.

### Stainless Steel Survival Kit for HPLC

The stainless steel survival kit contains a wide range of tubing, fittings, and tools necessary to set up and maintain your HPLC system: a selection of lengths and IDs of 1/16-inch tubing, plus nuts, ferrules, a ValvTool wrench, and a zero-dead-volume union.

Description	qty.	cat.#
Stainless Steel Survival Kit for HPLC	kit	25097

#### Contents:

Description	color	qty.	cat.#
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 5cm	red	3-pk.	25240
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 10cm	red	3-pk.	25241
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 20cm	red	3-pk.	25242
HPLC Capillary Tubing, SS, 1/16" x 0.005" x 30cm	red	3-pk.	25243
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 5cm	gray	3-pk.	25244
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 10cm	gray	3-pk.	25245
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 20cm	gray	3-pk.	25246
HPLC Capillary Tubing, SS, 1/16" x 0.007" x 30cm	gray	3-pk.	25247
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 5cm	blue	3-pk.	25248
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 10cm	blue	3-pk.	25249
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 20cm	blue	3-pk.	25250
HPLC Capillary Tubing, SS, 1/16" x 0.010" x 30cm	blue	3-pk.	25251
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 5cm	yellow	3-pk.	25252
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 10cm	yellow	3-pk.	25253
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 20cm	yellow	3-pk.	25254
HPLC Capillary Tubing, SS, 1/16" x 0.020" x 30cm	yellow	3-pk.	25255
1/16" Rheodyne® Style Nut		10-pk.	25095
1/16" Rheodyne® Style Ferrule		10-pk.	25096
ValvTool Wrench		ea.	25321
1/16" Stainless Steel Ferrules		10-pk.	20286
1/16" Stainless Steel Nuts		10-pk.	20287
Zero-Dead-Volume Internal Union		ea.	20148



Stainless Steel Survival Kit

### it's a fact

Survival Kits are convenient and economical, compared to ordering individual items!