



### Inert PEEK™ Tubing

- Replaces stainless steel, titanium, Teflon® or Tefzel® tubing.
- Less oxygen permeable and more temperature resistant (to 250°C) than Teflon® or Tefzel® tubing.
- Use with PEEK™ fingertight or flangeless fittings.
- Use to 7,000psi (48,263kPa).

Description	Color Code	qty.	cat.#
PEEK™ Tubing, 1/16" OD x 0.0025" ID x 1m	pink dash stripe	3-pk.	25320
PEEK™ Tubing, 1/16" OD x 0.005" ID x 3m	red stripe	ea.	25065
PEEK™ Tubing, 1/16" OD x 0.007" ID x 3m	yellow stripe	ea.	25066
PEEK™ Tubing, 1/16" OD x 0.010" ID x 3m	blue stripe	ea.	25067
PEEK™ Tubing, 1/16" OD x 0.020" ID x 3m	orange stripe	ea.	25068

### HPLC Stainless Steel Capillary Tubing

- 316 grade stainless steel.
- Precise pre-cut lengths.
- Smooth surface finish.
- Ultra clean.

Whether you need to replace system tubing as part of your troubleshooting or are looking to reduce the dwell volume of your system as you move to narrower columns, Restek has the quality tubing in the lengths and IDs you need. Each ID is color-coded so it is easy to identify and replace correctly.

Length	qty.	0.005" ID (red) cat.#	0.007" ID (gray) cat.#	0.010" ID (blue) cat.#	0.020" ID (yellow) cat.#
1/16" OD Tubing					
5cm	3-pk.	25240	25244	25248	25252
10cm	3-pk.	25241	25245	25249	25253
20cm	3-pk.	25242	25246	25250	25254
30cm	3-pk.	25243	25247	25251	25255

### Teflon® Tubing

- Ideal for mobile phase inlet lines.
- Chemically inert.
- Use to 500psi (3,447kPa) and 80°C.

Description	qty.	cat.#
Teflon® Tubing, 1/8" OD x 0.063" ID x 3m (1.6mm ID)	3m	25306
Teflon® Tubing, 1/8" OD x 0.094" ID x 3m (2.4mm ID)	3m	25307

### Tubing Clip

Securely holds 1/16-inch or 1/8-inch tubing in beaker, flask, or bottle up to 4mm thick.

Description	qty.	cat.#
Tubing Clip	5-pk.	25310