

Syringe Filters

- Variety of filter types, porosities, and diameters.
- Color coded for easy identification.
- Reusable storage container for 100-pk.
- 500- and 1,000- pack sizes for greater savings.



Size	Porosity	Color	100-pk. cat.#	500-pk. cat.#	1,000-pk. cat.#
Nylon					
13mm	0.22µm	yellow	26146		26200
13mm	0.45µm	pink	26147		26201
25mm	0.22µm	yellow	26148	26202	
25mm	0.45µm	pink	26149	26203	
PTFE (polytetrafluoroethylene)					
13mm	0.22µm	purple	26142		26196
13mm	0.45µm	orange	26143		26197
25mm	0.22µm	purple	26144	26198	
25mm	0.45µm	orange	26145	26199	
PVDF (polyvinylidene difluoride)					
13mm	0.22µm	black	26150	new!	26154
13mm	0.45µm	red	26151		
25mm	0.22µm	black	26152		
25mm	0.45µm	red	26153		

Nylon, PVDF - hydrophilic applications

PTFE - hydrophobic applications

free sample

Free sample packs of syringe filters are available by adding -248 to the catalog number.

Resprep™ SPE Tubes, Frits, & Caps

Empty Tubes (polypropylene)	Volume	qty.	cat.#
	1mL	50-pk.	26010
	3mL	50-pk.	26011
	6mL	50-pk.	26012
	15mL	50-pk.	26013
	sample reservoir, 25mL	12-pk.	26014
	sample reservoir, 70mL	12-pk.	26015
Frits (polyethylene)	Fits Tube Volume, Diameter	qty.	cat.#
	1mL, 6mm	100-pk.	26016
	3mL, 9mm	100-pk.	26017
	6mL, 1.2cm	100-pk.	26018
	15mL, 1.6cm	100-pk.	26019
	25mL, 2.0cm*	100-pk.	26020
	70mL, 2.6cm**	100-pk.	26021
Tube Caps	Fits Tube Volume	qty.	cat.#
	1mL	12-pk.	26001
	3mL	12-pk.	26002
	6mL	12-pk.	26003
	15mL	12-pk.	26004
	25mL*	12-pk.	26005
	70mL**	12-pk.	26006



Empty Tube



Frits



Tube Caps

Resprep™ SPE Connectors & End Caps

Connectors	Fits Tube Volume	qty.	cat.#
	1, 3, 6, 15mL	12-pk.	26007
	25mL*	12-pk.	26008
	70mL**	12-pk.	26009
Female Luer End Caps	universal	12-pk.	26000



Connectors

*For 20mL packed tubes.

**For 60mL packed tubes.