






Siltek®/Sulfinert® Treated and Silcosteel®-CR Treated Parker Fittings

A broad line of 1/16", 1/8" and 1/4" fittings is available with Siltek®/Sulfinert® or Silcosteel®-CR treatment. Because of expanding applications for these coatings, we have received many requests for a broader product offering. If you do not see everything you need, contact us for information on custom coating services.

Fitting Type	Size	Parker #	Siltek®/Sulfinert® Treated*		Silcosteel®-CR Treated	
			qty.	cat.#	qty.	cat.#
 Union	1/16"	1 SC 1	ea.	22520	ea.	22863
	1/8"	2 SC 2	ea.	22521	ea.	22864
	1/4"	4 SC 4	ea.	22522	ea.	22865
 Tee	1/16"	1 ET 1	ea.	22526	ea.	22866
	1/8"	2 ET 2	ea.	22527	ea.	22867
	1/4"	4 ET 4	ea.	22528	ea.	22868
 Reducing Union	1/8" to 1/16"	2 RU 1	ea.	22523	ea.	22869
	1/4" to 1/16"	4 RU 1	ea.	22524	ea.	22870
	1/4" to 1/8"	4 RU 2	ea.	22525	ea.	22871
 Elbow	1/16"	1 EE 1	—	—	ea.	22874
	1/8"	1 EE 2	ea.	22530	ea.	22875
	1/4"	1 EE 4	ea.	22531	ea.	22876
 Plug	1/16"	1 BLP 1	ea.	21539	ea.	22877
	1/8"	2 BLP 2	ea.	21540	ea.	22878
	1/4"	4 BLP 4	ea.	21541	ea.	22879
 Cross	1/8"	2 ECR 2	ea.	21542	ea.	22872
	1/4"	4 ECR 4	ea.	21543	ea.	22873
 Tube End Reducer	1/8" tube to 1/16"	2 TUR 1	ea.	21544	ea.	22880
	1/4" tube to 1/16"	4 TUR 1	ea.	21545	ea.	22881
	1/8" tube to 1/4"	2 TUR 4	ea.	21546	ea.	22882
	1/4" tube to 1/8"	4 TUR 2	ea.	21547	ea.	22883
 Port Connector	1/8"	2 PC 2	ea.	21548	ea.	22884
	1/4"	4 PC 4	ea.	21549	ea.	22885
	1/8" tube to 1/4"	2 PC 4	ea.	21550	ea.	22886
 Male Connector	1/8" to 1/8" NPT	2 MSC 2N	ea.	21551	ea.	22887
	1/4" to 1/4" NPT	4 MSC 4N	ea.	21552	ea.	22888
	1/16" to 1/8" NPT	1 MSC 2N	ea.	21553	ea.	22889
	1/8" to 1/4" NPT	2 MSC 4N	ea.	21554	ea.	22890
	1/4" to 1/8" NPT	4 MSC 2N	ea.	21555	ea.	22891
 Female Connector	1/8" to 1/8" NPT	2 FSC 2N	ea.	21556	ea.	22892
	1/4" to 1/4" NPT	4 FSC 4N	ea.	21557	ea.	22893
	1/4" to 1/8" NPT	4 FSC 2N	ea.	21558	ea.	22894
	1/8" to 1/4" NPT	2 FSC 4N	ea.	21559	ea.	22895
 Plug Valve, 2-Way	1/8"	2A PR4 VT SS	ea.	21586	—	—
	1/4"	4A PR4 VT SS	ea.	21587	—	—
 Ball Valve, 2-Way	1/8"	2A B2LJ2 SSP	ea.	21588	—	—
	1/4"	4A B2LJ2 SSP	ea.	21589	—	—

also available

For nontreated brass and stainless steel fittings, see pages 252-253.

new!

For the newest developments, visit our Performance Coatings Division website: www.restekcoatings.com

*See Frequently Asked Questions on page 391.

Please note: Nuts and ferrules are not treated unless requested (custom parts). Nuts and ferrules normally are not in contact with samples, and thus do not require coating.