

new!



method applications

Method	Application
US EPA	TO-17
ASTM	D-6196
NIOSH	2549
DIN EN ISO	16017

Specifications

Dimensions: 1/4" OD x 3-1/2" long
 Low sampling rates:
 0.01-0.20 L/min.
 (<10L total volume)
 Long-term storage caps are supplied with conditioned tubes

Thermal Desorption Unit (TDU) Tubes

- Variety of sorbents to collect a wide range of VOCs.
- Use glass tubes for maximum inertness in active sampling.
- Choose stainless steel tubes for either active or passive sampling. No sampling pump necessary for passive sampling with diffusion caps!
- Individually etched with unique serial number for convenient sample identification.
- Available unconditioned or preconditioned and ready to sample. Tubes are reusable after thermal desorption.

High-quality thermal desorption tubes by Markes International are now available from Restek. These sorbent tubes are suitable for ppt to ppm concentrations of volatile organic compounds (VOCs) in ambient, indoor, and industrial hygiene environments. Available in both stainless steel and glass (for thermally labile VOCs), they fit Markes ULTRA-UNITY, PerkinElmer, and Shimadzu thermal desorbers. Packed tubes come with a report detailing the total mass of sorbent in the tube; conditioned tubes also include a blank chromatogram.

Thermal Desorption Tube Sorbent	Applications
Tenax TA	Vapour phase organics from C6/7 to C26
Graphitized Carbon	Vapour phase organics from C5/6 to C14
Tenax GR/Carbopack™ B	Vapour phase organics from n-C5/6 to n-C20 (EPA Methods TO-14/TO-15/TO-17)
Carbopack™ B/Carbosieve™ SIII	Vapour phase organics from n-C2/3 to n-C12/14 (EPA Methods TO-14/TO-15/TO-17)
Tenax TA/Graphitized Carbon/Carboxen™ 1000	Vapour phase organics from C2/3 to C20
Carbopack™ C/Carbopack™ B/Carbosieve™ SIII	Vapour phase organics from n-C2/3 to n-C16/20 (EPA Methods TO-14/TO-15/TO-17)

Thermal Desorption Unit Tubes, Unconditioned and Conditioned & Capped

Description	qty.	Unconditioned		Conditioned & Capped	
		Stainless Steel	Glass	Stainless Steel	Glass
		cat.#	cat.#	cat.#	cat.#
TDU Tubes, Tenax TA	10-pk.	24056	24062	24080	24086
TDU Tubes, Graphitized Carbon	10-pk.	24057	24063	24081	24087
TDU Tubes, Tenax GR/Carbopack™ B	10-pk.	24058	24064	24082	24088
TDU Tubes, Carbopack™ B/Carbosieve™ SIII	10-pk.	24059	24065	24083	24089
TDU Tubes, Tenax TA/Graphitized Carbon/Carboxen™ 1000	10-pk.	24060	24066	24084	24090
TDU Tubes, Carbopack™ C/Carbosieve™ SIII	10-pk.	24061	24067	24085	24091



Stainless Steel, Conditioned and Capped



Glass, Unconditioned



Stainless Steel, Unconditioned

Thermal Desorption Unit Tubes, Empty

- Empty tubes for direct desorption of VOCs in liquids, solids, or pastes.
- Stainless steel: front sorbent retaining gauze fitted, rear gauze and gauze retaining spring supplied.
- Glass: with glass frit positioned 15mm from sampling end.

Description	qty.	Stainless Steel	Glass
		cat.#	cat.#
TDU Tubes, Empty	10-pk.	24054	24055

Thermal Desorption Unit Tubes, Calibration

Description	qty.	Stainless Steel	Glass
		cat.#	cat.#
TDU Tubes, Calibration, Tenax TA 1cm Bed	10-pk.	24075	24076
Description		qty.	cat.
Calibration Solution Loading Rig		ea.	24077
Calibration Solution Loading Rig 9.5mm Replacement Septa		10-pk.	24078
Certified Reference Standard, 100ng BTX on Tenax TA		10-pk.	24079

Thermal Desorption Unit Tubes, Accessories

Description	Benefits/Uses	qty.	cat.
1/4" Brass Cap and PTFE Ferrules	Use for long-term storage of blank/sampled tubes.	20-pk.	24068
1/4" PTFE Ferrules	Long-term storage caps.	20-pk.	24069
CapLok Tool	Use for tightening long-term storage caps.	ea.	24070
Pen Clip		10-pk.	24071
TubeMate Tool	Assists with tube packing.	ea.	24072
1/4" Stainless Steel Union and PTFE Ferrules	Use for connecting tubes in series.	10-pk.	24073
Diffusion Caps	Required for diffusive sampling with stainless steel tubes.	10-pk.	24074



Diffusion Caps