

Gas Sampling Bags

new!



Cali-5-Bond™ Gas Sampling Bags

- Totally nonpermeable and opaque, providing UV protection.
- Chemically inert—extremely rugged and portable.
- Extra strength—5 mil (0.14mm) thick.
- Easy to use.

Cali-5-Bond™ air and gas sampling bags provide a simple, reliable, and economic method of collecting air, gas, and liquid samples. The 5-layer construction (made by a patented process) ensures the physical integrity of any sample taken, providing a truly representative sample of the collection environment. Both grab and time-integrated samples can be taken with the use of a sampling pump. The twist-type valve with hose-barb connection enables secure attachment of $\frac{3}{16}$ " ID sample tubing. The septum port allows easy access via a gas-tight syringe. Bags should not be used at temperatures above 50°C (125°F) and should never be over inflated.



Description	qty.	cat.#
0.5L 6" x 8"	5-pk.	24092
1L 8" x 8"	5-pk.	24093
2L 8" x 12"	5-pk.	24094
5L 8" x 23"	5-pk.	24095
10L 16" x 15"	5-pk.	24096
22L 16" x 25"	5-pk.	24097
44L 24" x 25"	5-pk.	24098

Tedlar® Sampling Bags

- Find the bags you need—we offer sizes from 0.5 liters to 100 liters.
- Unique all-in-one septum and valve fitting make these lightweight and easy to use.
- Polypropylene or stainless steel valve.

The unique design of these Tedlar® sample bags incorporates the sampling septum directly in the valve, providing easier use and lighter weight than other styles. We offer two types of bags: one with a polypropylene valve and one with a stainless steel valve. Both valves conveniently connect to $\frac{3}{16}$ " ID Teflon® tubing.

Description	qty.	Polypropylene Valve	Stainless Steel Valve
		cat.#	cat.#
0.5L 6" x 6"	10-pk.	22049	22038
1L 7" x 7"	10-pk.	22050	22039
3L 9.5" x 10"	10-pk.	22051	22040
5L 12" x 12.5"	10-pk.	22052	22041
10L 11.75" x 22"	10-pk.	22053	22042
12L 13" x 24"	10-pk.	22054	22043
25L 17.5" x 24"	5-pk.	22055	22044
40L 24" x 24.25"	5-pk.	22056	22045
80L 28.25" x 32.5"	5-pk.	22057	22046
100L 28" x 36"	3-pk.	22058	22047
Replacement Septum	10-pk.	22059	22048