

Natural Gas and Refinery Gas Standards

Natural Gas and Refinery Gas Standards

- Each available in three varying concentrations.
- Mini-regulator designed specially for these standards.



Silvia Martinez
Innovations Chemist
5+ years of service!

Natural Gas Standards

Available in three mixes, from lean to rich. Each has an extended list of C6+ components.

	Natural Gas Standard #1 cat.# 34438, ea. % each compound**	Natural Gas Standard #2 cat.# 34439, ea. % each compound**	Natural Gas Standard #3 cat.# 34440, ea. % each compound**
nitrogen	1.000	2.500	5.000
carbon dioxide	0.500	1.000	1.500
methane UHP	94.750	85.250	70.000
ethane UHP	2.000	5.000	9.000
propane	0.750	3.000	6.000
isobutane	0.300	1.000	3.000
n-butane	0.300	1.000	3.000
isopentane	0.150	0.500	1.000
n-pentane	0.150	0.500	1.000
hexanes plus EX2*	0.100	0.250	0.500
Concentration	mole	mole	mole
Volume	13.16L @ 200psig	13.16L @ 200psig	5.5L @ 75psig
Ideal Heating Value (Dry BTU/SCF)	1048 gross	1142 gross	1317 gross

*Contact Restek or your Restek representative for a complete list of hexanes plus EX2.

**Precise concentrations are provided on the data sheet included with each cylinder and may vary slightly from those listed here.

Refinery Gas Standards

Available in three mixes with varying C5 unsaturates or extended C6+ components.

	Refinery Gas Standard #1 cat.# 34441, ea. % each compound**	Refinery Gas Standard #2 cat.# 34442, ea. % each compound**	Refinery Gas Standard #5 cat.# 34443, ea. % each compound**
hydrogen	40.750	12.500	12.500
argon	0.500	1.000	1.000
nitrogen	4.000	37.200	37.200
carbon monoxide	1.000	1.000	1.000
carbon dioxide	3.000	3.000	3.000
methane	8.500	5.000	5.000
ethane	6.000	4.000	4.000
ethylene	2.000	2.000	2.000
acetylene	—	1.000	1.000
propane	7.000	6.000	6.000
propylene	3.000	3.000	3.000
propadiene	0.850	1.000	1.000
cyclopropane	—	0.040	—
isobutane	6.000	5.000	5.000
n-butane	4.000	4.000	4.000
isobutylene	2.000	1.000	1.000
1,3 butadiene	3.000	3.000	3.000
cis-2-butene	2.000	2.000	2.000
trans-2-butene	2.000	3.000	3.000
butene-1	2.000	2.000	2.000
2-methyl-2-butene	—	0.200	0.200
isopentane	1.000	1.000	1.000
n-pentane	1.000	1.000	1.000
cis-2-pentene	—	0.400	0.400
trans-2-pentene	—	0.160	0.200
pentene-1	—	0.400	0.400
n-hexane	0.500	0.100	—
hexanes plus EX	—	—	0.100
Concentration	mole	mole	mole
Volume	5.2L @ 70psig	4.9L @ 60psig	4.6L @ 60psig

**Precise concentrations are provided on the data sheet included with each cylinder and may vary slightly from those listed here.

Mini-Regulator for natural gas and refinery gas standards

- 0–300psig inlet pressure range.
- 0–15psig outlet pressure range.
- Supplied with 0–15psig outlet pressure gauge, brass CGA 170 nut and nipple.

Description
Mini-Regulator

qty.
ea.

cat.#
22032



cylinder
design

DCG Partnership Cylinders:

Size: 7.6 x 24 cm.

CGA-170/110 connection.

U.S. D.O.T. Specs:

DOT-4B-240ET

Please note: This cylinder is not approved for use in Canada.

