

Method 529 (Nitroaromatics & Nitramines)

Nitroaromatics and Nitramine Explosives in Drinking Water (14 components)

3,5-dinitroaniline	2-nitrotoluene
1,3-dinitrobenzene	3-nitrotoluene
2-amino-4,6-dinitrotoluene	4-nitrotoluene
4-amino-2,6-dinitrotoluene	RDX
2,4-dinitrotoluene	tetryl
2,6-dinitrotoluene	1,3,5-trinitrobenzene
nitrobenzene	2,4,6-trinitrotoluene

1,000µg/mL each in acetonitrile, 1mL/ampul
cat. # 33900 (ea.)

3,4-Dinitrotoluene

3,4-dinitrotoluene
2,000µg/mL in ethyl acetate, 1mL/ampul
cat. # 33901 (ea.)

2-Nitromesitylene

2-nitromesitylene
2,000µg/mL in methanol, 1mL/ampul
cat. # 33902 (ea.)

1,2,4-Trimethyl-5-nitrobenzene

1,2,4-trimethyl-5-nitrobenzene
2,000µg/mL in methanol, 1mL/ampul
cat. # 33903 (ea.)

Nitrobenzene-d5

nitrobenzene-d5
2,000µg/mL in methylene chloride, 1mL/ampul
cat. # 33904 (ea.)

Method 531.1, 531.2 (Carbamates)

Internal Standard

4-bromo-3,5-dimethylphenyl-N-methylcarbamate (BDMC)
100µg/mL in methanol, 1mL/ampul
cat. # 32274 (ea.)

531.1 Performance Check Mix

aldicarb sulfoxide	100µg/mL	3-hydroxycarbofuran	2
BDMC	10	methiocarb	20

In methanol, 1mL/ampul
cat. # 32275 (ea.)

531.1 Carbamate Pesticide Calibration Mixture

(10 components)	
aldicarb	3-hydroxycarbofuran
aldicarb sulfone	methiocarb
aldicarb sulfoxide	methomyl
carbaryl (Sevin®)	oxamyl
carbofuran	propoxur (Baygon®)

100µg/mL each in methanol, 1mL/ampul
cat. # 32273 (ea.)

531.2 Carbamate Pesticide Calibration Mixture

(11 components)	
aldicarb	methiocarb
aldicarb sulfone	methomyl
aldicarb sulfoxide	1-naphthol
carbaryl (Sevin®)	oxamyl
carbofuran	propoxur (Baygon®)
3-hydroxycarbofuran	

100µg/mL in acetonitrile, 1mL/ampul
cat. # 32435 (ea.)

Method 532 (Phenylurea Pesticides)

Phenylurea Surrogate Mixture

carbazole	monuron
-----------	---------

500µg/mL each in methanol:acetonitrile (50:50), 1mL/ampul
cat. # 32433 (ea.)

Phenylurea Pesticide Mixture (8 components)

diflubenzuron	propanil
diuron	siduron
fluometuron	tebuthiuron
linuron	thidiazuron

200µg/mL each in acetonitrile:acetone (90:10), 1mL/ampul
cat. # 32434 (ea.)

also available

See the LC Applications section for a carbamate pesticides chromatogram - [page 532](#).

did you know?

We have more than 2,000 pure, characterized, neat compounds in our inventory! If you do not see the EXACT mixture you need listed on any of these pages, call us.

See [page 427](#) for our Custom Reference Materials Request Form.