

Pennsylvania

PA DEP UST Standard (11 components)

benzene	naphthalene
1,2-dibromoethane	toluene
1,2-dichloroethane	<i>m</i> -xylene
ethylbenzene	<i>o</i> -xylene
isopropyl benzene	<i>p</i> -xylene
methyl <i>tert</i> -butyl ether (MTBE)	

2,000µg/mL each in P&T methanol, 1mL/ampul
cat. # 30433 (ea.)

Tennessee/Mississippi

DRO Mix (Tennessee/Mississippi) (16 components)

<i>n</i> -decane (C10)	<i>n</i> -octadecane (C18)
<i>n</i> -undecane (C11)	<i>n</i> -nonadecane (C19)
<i>n</i> -dodecane (C12)	<i>n</i> -eicosane (C20)
<i>n</i> -tridecane (C13)	<i>n</i> -heneicosane (C21)
<i>n</i> -tetradecane (C14)	<i>n</i> -docosane (C22)
<i>n</i> -pentadecane (C15)	<i>n</i> -tricosane (C23)
<i>n</i> -hexadecane (C16)	<i>n</i> -tetracosane (C24)
<i>n</i> -heptadecane (C17)	<i>n</i> -pentacosane (C25)

1,000µg/mL each in methylene chloride, 1mL/ampul
cat. # 31214 (ea.)

Gasoline Component Standard (10 components)

benzene	500µg/mL	1,2,4-trimethylbenzene	1000
ethylbenzene	500	2,2,4-trimethylpentane	1500
heptane	500	<i>m</i> -xylene	1000
2-methylpentane	1500	<i>o</i> -xylene	1000
toluene	1500	<i>p</i> -xylene	1000

10,000µg/mL total in P&T methanol, 1mL/ampul
cat. # 30486 (ea.)



Don't see the UST mix you need? We can custom blend a UST mix to meet the requirements of your method. Use the form on [page 427](#) or visit our web site at www.restek.com.

Texas

Texas TNRCC Method 1006

TNRCC 1006 Retention Time Marker Mix

(9 components)

<i>n</i> -hexane (C6)	<i>n</i> -hexadecane (C16)
<i>n</i> -heptane (C7)	<i>n</i> -heneicosane (C21)
<i>n</i> -octane (C8)	<i>n</i> -octacosane (C28)
<i>n</i> -decane (C10)	<i>n</i> -pentatriacontane (C35)
<i>n</i> -dodecane (C12)	

200µg/mL in pentane, 1mL/ampul
cat. # 31814 (ea.)

Texas TNRCC Method 1005

TNRCC 1005 Retention Time Markers Mix

(9 components)

<i>n</i> -hexane (C6)	<i>n</i> -octacosane (C28)
<i>n</i> -dodecane (C12)	<i>n</i> -pentatriacontane (C35)

200µg/mL each in pentane, 1mL/ampul
cat. # 31698 (ea.)

TX TPH Locator Mix

<i>n</i> -hexane (C6)	<i>n</i> -octacosane (C28)
<i>n</i> -decane (C10)	

200µg/mL each in pentane, 1mL/ampul
cat. # 31482 (ea.)

TX TPH Calibration Mix

diesel fuel #2 composite	unleaded gasoline composite
--------------------------	-----------------------------

10,000µg/mL each in pentane, 1mL/ampul
cat. # 31483 (ea.)

TX TPH Matrix Spike Mix

diesel fuel #2 composite	unleaded gasoline composite
--------------------------	-----------------------------

10,000µg/mL each in P&T methanol, 1mL/ampul
cat. # 31484 (ea.)

Alternate Boiling Point/Carbon Number

Distribution Marker Stock Standard (9 components)

<i>n</i> -hexane (C6)	<i>n</i> -heneicosane (C21)
<i>n</i> -octane (C8)	<i>n</i> -octacosane (C28)
<i>n</i> -decane (C10)	<i>n</i> -pentatriacontane (C35)
<i>n</i> -dodecane (C12)	<i>n</i> -hexatriacontane (C36)
<i>n</i> -hexadecane (C16)	

200µg/mL each in pentane, 1mL/ampul
cat. # 31639 (ea.)

 α, α, α -Trifluorotoluene

2,000µg/mL in P&T methanol, 1mL/ampul	cat. # 30048 (ea.)
2,500µg/mL in P&T methanol, 1mL/ampul	cat. # 30068 (ea.)
10,000µg/mL in P&T methanol, 1mL/ampul	cat. # 30083 (ea.)

1-Chlorooctane

10,000µg/mL in P&T methanol, 1mL/ampul	cat. # 30084 (ea.)
--	--------------------

1-Chlorooctadecane Mix

1-chlorooctadecane	
10,000µg/mL in methylene chloride, 1mL/ampul	cat. # 31098 (ea.)

free literature

Texas UST Monitoring

Download your free copy from www.restek.com.

Fast Facts

lit. cat. # 59394

did you know?

We have more than 2,000 pure, characterized, neat compounds in our inventory! If you do not see the EXACT mixture you need listed on any of these pages, call us.

See [page 427](#) for our Custom Reference Materials Request Form.