

also available

Custom total sulfur & total nitrogen in isooctane check standards also available.

Call our Analytical Reference Department at 800-356-1688 or 814-353-1300, or your Restek representative for details.

free data

Available on Our Website: Lot Certificates, Data Packs, and MSDSs

For complete information detailing manufacturing and testing for Restek inventoried reference standards, visit our website at www.restek.com.

To view lot certificates and/or an MSDS, enter the catalog number of the product in the Search feature. For a free data pack (Adobe® PDF file), enter the catalog number and lot number of the product.

Sulfur in Isooctane Calibration Kits and Check Standards *cont'd*

Cal Kit SISO 2.5 - 50 ppm

kit

Blank

2.50 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
5.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
10.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
15.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
20.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
25.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
50.00 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane

Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane 2.5-50 ppm range. Set of eight 1mL pre-scored ampuls.
cat. # 33037 (kit)

Check Standard SISO 30

30 ppm total sulfur by weight from di-*n*-butylsulfide in isooctane. Set of five 1mL pre-scored ampuls.

cat. # 33038 (ea.)

Cal Kit SISO 50 - 1000 ppm

kit

Blank

50 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
75 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
100 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
250 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
500 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
1000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane

Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane 50-1000 ppm range. Set of seven 1mL pre-scored ampuls.
cat. # 33039 (kit)

Check Standard SISO 300

300 ppm total sulfur by weight from di-*n*-butylsulfide in isooctane. Set of five 1mL pre-scored ampuls.

cat. # 33040 (ea.)

Cal Kit SISO 1000 - 6000

kit

Blank

1000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
1500 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
2000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
4000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane
6000 w/w ppm total sulfur from di-*n*-butylsulfide in isooctane

Calibration kit for total sulfur by weight from di-*n*-butylsulfide in isooctane 1000 - 6000 ppm range. Set of six 1mL pre-scored ampuls.
cat. # 33041 (kit)

Check Standard SISO 3000

3000 ppm total sulfur by weight from di-*n*-butylsulfide in isooctane. Set of five 1mL pre-scored ampuls.

cat. # 33042 (ea.)

Total Sulfur & Total Nitrogen in Isooctane Calibration Kits **new!**

ASTM Methods D3120, D4045, D4629, D5453, D5762, D6069, D6920

Cal Kit SNISO 0.125 - 5.0

kit

Blank

0.125 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
0.25 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
0.50 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
1.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
2.50 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
5.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of seven 1mL pre-scored amber ampuls.
cat. # 33057 (kit)

Cal Kit SNISO 5.0 - 50.0

kit

Blank

5.00 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
10.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
25.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
50.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of five 1mL pre-scored amber ampuls.
cat. # 33058 (kit)

Cal Kit SNISO 50.0 - 1000

kit

Blank

50.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
75.0 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
100 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
250 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
500 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane
1000 w/w ppm total sulfur from thiophene & total nitrogen from pyridine in isooctane

Set of seven 1mL pre-scored amber ampuls.
cat. # 33059 (kit)