

USP <467>

The United States Pharmacopeia (USP) general chapter <467> *Residual Solvents* is a widely used compendial method intended for identifying and quantifying residual solvents in drug substances, drug products, and excipients. In an attempt to better mirror the International Conference on Harmonization (ICH) guidelines, the USP is proposing a more comprehensive methodology in residual solvent testing—the current USP30/NF25. The ICH publishes a guideline (Q3C) listing the acceptable amounts of solvent residues that can be present. In the ICH guideline, residual solvents are summarized by class, according to their toxicity. Class 1 compounds are carcinogenic compounds that pose a risk to both the consumer and the environment. The use of these solvents is to be avoided, but if they are used, they must be tightly controlled. Class 2 compounds are nongenotoxic animal carcinogens, and concentrations of these compounds should be limited. Chromatographic analysis is needed for both the Class 1 and Class 2 residual solvents.

Initially set to become effective July 1, 2007, the implementation of the current version of USP <467>: *Residual Solvents* has been delayed until July 1, 2008. Until that time, the *Other Analytical Procedures* section of the previous version will be retained.

also available

For a list of OVI retention times, see page 721.

USP <467> Singles

Volume is 1mL/ampul.

Compound	Solvent	Conc.	cat.# (ea.)
acetonitrile	DMSO	2.05mg/mL	36281
benzene	DMSO	10mg/mL	36282
carbon tetrachloride	DMSO	20mg/mL	36283
chlorobenzene	DMSO	1.8mg/mL	36284
chloroform	DMSO	0.3mg/mL	36285
cyclohexane	DMSO	19.4mg/mL	36286
1,1-dichloroethene	DMSO	40mg/mL	36287
1,2-dichloroethane	DMSO	25mg/mL	36288
cis-1,2-dichloroethylene	DMSO	4.67mg/mL	36289
trans-1,2-dichloroethylene	DMSO	4.67mg/mL	36290
1,2-dimethoxyethane	DMSO	0.5mg/mL	36291
N,N-dimethylacetamide	DMSO	5.45mg/mL	36292
N,N-dimethylformamide	DMSO	4.4mg/mL	36293
1,4-dioxane	DMSO	1.9mg/mL	36294
2-ethoxyethanol	DMSO	0.8mg/mL	36295
ethylbenzene	DMSO	1.84mg/mL	36296
ethylene glycol	DMSO	3.1mg/mL	36297
formamide	DMSO	1.1mg/mL	36298
hexane	DMSO	1.45mg/mL	36299
methanol	DMSO	15mg/mL	36401
2-methoxyethanol	DMSO	0.25mg/mL	36402
methylbutylketone	DMSO	0.25mg/mL	36400
methylcyclohexane	DMSO	5.9mg/mL	36403
methylene chloride	DMSO	3mg/mL	36404
N-methylpyrrolidone	DMSO	2.65mg/mL	36405
nitromethane	DMSO	0.25mg/mL	36406
pyridine	DMSO	1mg/mL	36407
sulfolane	DMSO	0.8mg/mL	36413
tetrahydrofuran	DMSO	3.6mg/mL	36408
tetralin	DMSO	0.5mg/mL	36409
toluene	DMSO	4.45mg/mL	36410
1,1,1-trichloroethane	DMSO	50mg/mL	36411
trichloroethene	DMSO	0.4mg/mL	36412
m-xylene	DMSO	6.51mg/mL	36414
o-xylene	DMSO	0.97mg/mL	36415
p-xylene	DMSO	1.52mg/mL	36416

DMSO = dimethyl sulfoxide

These mixtures represent the changes made in USP30/NF25 effective July 1, 2008.

Residual Solvents - Class 1

new!

benzene	10mg/mL	1,1-dichloroethene	40
carbon tetrachloride	20	1,1,1-trichloroethene	50
1,2-dichloroethane	25		

In dimethyl sulfoxide, 1mL/ampul
cat. # 36279 (ea.)

These mixtures represent the changes made in USP30/NF25 effective July 1, 2008. *cont'd*

Residual Solvents Class 2 - Mix B

new!

(8 components)

chloroform	60µg/mL	nitromethane	50
1,2-dimethoxyethane	100	pyridine	200
n-hexane (C6)	290	tetralin	100
2-hexanone	50	trichloroethene	80

In dimethyl sulfoxide, 1mL/ampul
cat. # 36280 (ea.)

also available

For other reference mixes for USP <467> and European Pharmacopoeia, see page 510.

These mixtures represent the changes made in USP28/NF23 effective January 1, 2005.

Residual Solvents Class 2 - Mix A

(15 components)

acetonitrile	2.05mg/mL	methylcyclohexane	5.90
chlorobenzene	1.80	methylene chloride	3.00
cyclohexane	19.40	tetrahydrofuran	3.45
cis-1,2-dichloroethene	4.70	toluene	4.45
trans-1,2-dichloroethene	4.70	m-xylene	6.51
1,4-dioxane	1.90	o-xylene	0.98
ethylbenzene	1.84	p-xylene	1.52
methanol	15.00		

In dimethyl sulfoxide, 1mL/ampul
cat. # 36271 (ea.)

Class III solvents are available as custom mixes. See page 427 for our custom reference material request form.

Residual Solvents Class 2 - Mix B

(8 components)

chloroform	300µg/mL	nitromethane	250
1,2-dimethoxyethane	500	pyridine	1,000
n-hexane (C6)	1,450	tetralin	500
2-hexanone	250	trichloroethene	400

In dimethyl sulfoxide, 1mL/ampul
cat. # 36272 (ea.)

free literature

USP <467> Residual Solvents Method

Download your free copy from our website!

Flyer
lit. cat. # 580166

Residual Solvents Class 2 - Mix C

(8 components)

2-ethoxyethanol	800µg/mL	2-methoxyethanol (methyl Cellosolve®)	250
ethylene glycol	3,100	N-methylpyrrolidone	2,650
formamide	1,100	sulfolane	800
N,N-dimethylacetamide	5,450		
N,N-dimethylformamide	4,400		

In dimethyl sulfoxide, 1mL/ampul
cat. # 36273 (ea.)