

HPLC Applications by Column Phase (Alphabetical Order) - see pages 734-735 for our Applications by Compound Class Index

Allure® AK	
carbonyls	.321, 528
Allure® Basix	
acetaminophen, narcotic analgesics	.544
acetaminophen, pseudoephedrine	
& caffeine	.559
albuterol	.553
antiarrhythmics	.546, 547
antibiotics	.550
antidepressants	.550
herbicides	.531
ibuprofen & pseudoephedrine	.559
sedatives	.558
steroids	.556
Allure® Biphenyl	
antidepressants	.551
cannabinoids	.526
corticosteroids	.553, 554
NSAIDs	.557
steroids	.554, 555
Allure® C18	
antibiotics	.550
herbicides	.532
pesticides	.533
Allure® Organic Acids	
carboxylic acids	.540
fruit juice acids	.511, 540, 541
Allure® PFP Propyl	
antibiotics	.549
catecholamines	.525
cocaine, ecgonine methyl ester	.527
nucleic acid bases	.523
opiates & metabolites	.526
pesticides	.533
pHidelity® C18	
basic drugs	.556
serotonin reuptake inhibitors	.552
Pinnacle™ DB Biphenyl	
steroids	.555
Pinnacle™ DB C18	
parabens	.537
vanilla bean extracts	.541
xanthines	.558
Pinnacle™ DB PFP Propyl	
benzodiazepines	.527
Pinnacle™ II Amino	
lactulose concentrate	.558
maple syrup	.541
sugars	.541
Pinnacle™ II Biphenyl	
explosives	.317, 529
Pinnacle™ II C18	
allicin	.536
capsaicinoids	.539
morphine sulfate	.545
phenolic antioxidants	.536, 539
Pinnacle™ II C8	
parabens	.537
Pinnacle™ II Cyano	
corticosteroids	.553
piperine	.539
Pinnacle™ II PAH	
PAHs	.535
Pinnacle™ II Phenyl	
parabens	.537
sorbic & benzoic acids	.539
Pinnacle™ II Silica	
hydrocodone bitartrate	.545
tocopherols	.536
Ultra Aqueous C18	
amino acids (aromatics)	.524
analgesics	.544
carboxylic acids	.540
herbicides	.530, 531
phenethyl glucosinolate	.536
sudan dyes	.542
Ultra C18	
acetaminophen, narcotic analgesics	.544
aldehydes, ketones	.528
beclomethasone	.548
corticosteroids	.553
explosives	.529
herbicides	.530
hydrocodone bitartrate, acetaminophen	.559
nitrofurans metabolites	.543
vitamins (fat soluble)	.538
Ultra C8	
oxycodone	.545
vanillin & ethyl vanillin	.542
Ultra Carbamate	
pesticides (carbamates)	.532
Ultra Cyano	
antiarrhythmics	.545
Ultra IBD	
antiarrhythmics	.546
antibiotics	.548, 550
antidepressants	.551
cimetidine	.552
diflubenzuron (pesticide)	.530
glyburide	.552
nucleosides, nucleotides,	
& nucleic acid bases	.523
vitamins (water soluble)	.538
Ultra PFP	
digitalis extract	.546
levothyroxine	.556
purines, pyrimidines	.524
Ultra Phenyl	
guaifenesin, codeine	.559
Ultra Quat	
paraquat, diquat	.328, 534
Viva C18	
peptides	.524
proteins	.525

HPLC & GC Applications
Compound Index

starting on page 736

Searching for a chromatogram?
www.restek.com