

Custom Packed Columns

To order, specify the following:

- 1) column dimensions (length, ID) and tubing material
- 2) packing description (percent coating and phase, support mesh size, and treatment)
- 3) column configuration (instrument manufacturer, model number, on-column injection or not) and with or without nuts and ferrules

Ordering Example: (6' x 1/8") (stainless steel) (3%) (Rtx®-1) (Silcoport® 80/100) (Agilent 6890) (on-column injection) (fittings kit).

Please use the custom order form on page 142 or visit www.restek.com/packed



Custom Micropacked Columns

To order, contact your Restek representative and specify the following:

- 1) physical dimensions (length, OD, ID, and tubing material)
- 2) packing description (percent coating and phase, support mesh size)
- 3) installation kit (see page 133), frit type

Ordering Example: (2 m x 1/16" OD x 1.00 mm ID) (Siltek®-treated tubing) (5%) (Carbowax® 20M) (CarboBlack B) (80/120) (installation kit for valve applications, cat. #21065) (Siltek® frits)

Please use the custom order form on page 142 or visit www.restek.com/packed

did you know?

Packing material in packed and micropacked columns is secured using wire braids or frits. This prevents packing material from exiting the column.

Frits—A new alternative to glass wool and braided end plugs!

Hastelloy® and Siltek® treated frits are now available for select packed and micropacked columns!



Hastelloy® frit Siltek® frit

Fill out the form on page 142, visit www.restek.com/packed, contact Customer Service or your Restek representative for pricing and availability.

ordering note

For international pricing on custom packed or micropacked columns, please contact your Restek representative.